

# Particle Count Strategies on Alarm Settings

by Jason Kelly

April 2017

## Setting APPROPRIATE Alarm Limits in Sterile Manufacturing Processes.

Jason Kelly – Lighthouse Worldwide Solutions

### Introduction

Setting inappropriate Alert and Action alarms on your Environmental Monitoring System (EMS) could be a big contributor to loss revenue, longer downtimes and a lot of discarded product. Misinterpretation of cGMP guidelines and an inadequate approach to the use of your EMS alarming functionality may set you up for failure right from the start.

By partnering with the right EMS supplier you can have a major impact on your manufacturing process, product revenue and regulatory approval process. Installation and proper use of an EMS depends on investing in the right supplier and on them having adequate resources to assist in maintaining the system and providing the right technical support required to enable business continuity. PICs GMP Annex 11 Computerized Systems discusses the need for a strong vendor partnership and ongoing support. So selecting the right supplier will assist in a monitoring system working for you and not against you.

When setting EMS alarm limits in sterile manufacturing, companies are restricting themselves by following a table of recommended alarm limits based on a volume of air sampled per cubic meter adopted from ISO 14644-1 a standard with an emphasis on Cleanroom Certification and not on continuous particle monitoring. (Continuous particle monitoring means continuous during the production process time and not 24/7).

Why would you implement alert and action alarms based on a cleanroom certification process which is a snap shot of a set volume which correlates to a table of maximum permitted particles per m<sup>3</sup>. This quickly becomes an apple to oranges scenario when the medium to collect real-time particle counts is a remote particle counter with a sample rate based on a 60 second interval or 1 cubic foot per minute sample rate and the time frame is continuous for the duration of the manufacturing process. In short it does not make sense and is a contributing factor to nuisance alarms being triggered during the production process. PICs GMP also tabulates alarm limits for Grade A-D environments for classification which are also followed for real time monitoring.

GMP guidelines recommend data for sizes of 0.5µm and 5.0µm is collected. The 5.0µm size range being more stringent as that size and above is normally considered the size range of bacterial types (viable particles) of particles. Below are examples of the PICs GMP Annex 1 table and ISO-14644-1 tables with alarm limits based on a cubic meter sample volume. These tables are generally used to set alarm limits. You will notice the 5µm limit is 20 counts per sample volume of 1m<sup>3</sup> for PICs GMP and this was previously 29 counts in the ISO 14644-1 (1999) table and is now removed based on statistical limitations when sampling in low concentration areas such as Grade A or ISO Class 5 environments in ISO 14644-1 2015 revision.

**Note** with the low levels of 5µm particles normally observed in these environments statistically it is not achievable to get meaningful data to make appropriate classification or real time trends.

Grade	Maximum permitted number of particles per m <sup>3</sup> equal to or greater than the tabulated size			
	At rest		In operation	
	0.5 μm	5.0 μm	0.5 μm	5.0 μm
A	3 520	20	3 520	20
B	3 520	29	352 000	2 900
C	352 000	2 900	3 520 000	29 000
D	3 520 000	29 000	Not defined	Not defined

Fig 1 PICs GMP table from Annex 1 2008

ISO Classification Number(N)	0.1 μm	0.2 μm	0.3 μm	0.5 μm	1.0 μm	5.0 μm
ISO 1	b 10	d -2	d	d	d	e
ISO 2	100	24	10	d -4	d	e
ISO 3	1,000	237	102	35	d 8	e
ISO 4	10,000	2,370	1,020	352	83	e
ISO 5	100,000	23,700	10,200	3,520	832	d,e,f 29
ISO 6	1,000,000	237,000	102,000	35,200	8,320	293
ISO 7	c	c	c	352,000	83,200	2,930
ISO 8	c	c	c	3,520,000	832,000	29,300
ISO 9	c	c	c	35,200,000	8,320,000	293,000

5.0 μm should be zero according to ISO/DIS 14644-1.2 Table 1 notes;

d) Sampling and statistical limitations for particles in low concentrations make classification inappropriate

e) ... Greater than 1 micron particles make classification at this particle size inappropriate due to potential particle losses in sampling system

f) Specify particle size in association with ISO Class 5, the marcoparticle descriptor M may be adapted.

Fig 2 ISO 14644-1 TABLE (2015)

**Table 1 — Selected airborne particulate cleanliness classes for cleanrooms and clean zones**

ISO classification number (N)	Maximum concentration limits (particles/m <sup>3</sup> of air) for particles equal to and larger than the considered sizes shown below (concentration limits are calculated in accordance with equation (1) in 3.2)					
	0,1 μm	0,2 μm	0,3 μm	0,5 μm	1 μm	5 μm
ISO Class 1	10	2				
ISO Class 2	100	24	10	4		
ISO Class 3	1 000	237	102	35	8	
ISO Class 4	10 000	2 370	1 020	352	83	
ISO Class 5	100 000	23 700	10 200	3 520	832	29
ISO Class 6	1 000 000	237 000	102 000	35 200	8 320	293
ISO Class 7				352 000	83 200	2 930
ISO Class 8				3 520 000	832 000	29 300
ISO Class 9				35 200 000	8 320 000	293 000

NOTE Uncertainties related to the measurement process require that concentration data with no more than three significant figures be used in determining the classification level

Fig.3 ISO 14644-1 TABLE (2009)

For example following PICs GMP the limit set in a Grade A environment at  $0.5\mu\text{m}$  is  $3520\text{ particles}/\text{m}^3$  and  $5.0\mu\text{m}$  is  $20\text{ particles}/\text{m}^3$  and if using a “rolling average or summation” to track this alarm limit with  $5\mu\text{m}$  at 20 particles they are pretty tight limits and the way micro and production managers set the action alarms meant that if 1 sample period was  $>20$  particles then it was seen an “Action Alarm” referred to as an excursion.

An SOP is typically followed and that usually meant having to segregate the batch running through the filling machine and conduct a root cause investigation and then go into the cleanroom bringing in more people and equipment inadvertently more contamination. When the root cause analysis was completed 99% of the time no plausible root cause was found! However this process managed to add more contamination to the cleanroom environment and increase down time, driving up production costs. Sounds like a futile exercise, waste of resources chasing a non-event simply because of inappropriate alarm settings in the first place.

The above introduction helps build up the picture for the solutions outlined below. By the way those root cause analysis may cost upward from \$10-\$20K to investigate per event with labor, other resources and laboratory time costs not to mention sometimes discarded batches which could be a couple of \$100K. So why did the root cause fail to identify a problem? It’s simple - alarm limits were not set up based on the actual process and the baseline of that process and the risk to product contamination. But were set up based on a table for certification purposes. It was assumed that an action alarm meant that there was a failure of the clean air system, in other words it was not set up to look at adverse trends that may impact product quality in the actual process environment.

The problem is also that the user of the particle counter may not have correctly understood particle counter limitations and functions. Particle counters count all particles regardless of viable or non-viable particles and the sample probe needs to be positioned in a validated position to report meaningful data from the process. The main function of the particle counter is to alarm if the clean air flowing over the process is no longer acting as a clean sterile barrier based on trending data and repeated alarm events. In the case of sterile injectable or compounding applications, appropriate thought and validation needs to be undertaken to verify the best position of the sample probe, the set-points, SPC triggers for alert and action limits as well as adequate SOPs to be executed when trending alarms are activated.



Fig.4 Study using a portable particle counter to validate sample probe placement and height inside LAF cabinet based on process manipulations

## Sample Probe Location, Risk Assessment & Setting Alarm limits

### (1) Probe Placement

The sample probe needs to be positioned above the process zone typically within a foot. (FDA guidelines). To correctly identify the sample position in a BSC Cabinet or filling line the actual physical activities in those environments need to be considered as well as the process being monitored.




Fig 5 & 6 Example of a remote particle counter sample probe placement on filling line as sterile vials exit a sterilization oven. The remote particle counter is under the filling machine and the tubing is kept to a minimum without any bends if possible.

The sample probe should face towards the incoming clean air from the HEPA filter in that environment. A stainless steel probe is the best type of probe material to use. It should also be lined internally with Bev-A-Line® tubing which has excellent particle anti-static and transportation properties. Care must be taking to limit the length of the tubing and any bends introduced. Placing the remote particle counter as close to the clean zone to be monitored is the best approach and again stainless steel is the preferred option for remote particle counter housing. PICs GMP Annex 1 recommends keeping sample tubing as short as possible. Section 11 below confirms the importance of sample tubing and sample probe placement.

***11. Airborne particle monitoring systems may consist of independent particle counters; a network of sequentially accessed sampling points connected by a manifold to a single particle counter; or a combination of the two. The system selected must be appropriate for the particle size considered. Where remote sampling systems are used, the length of tubing and the radii of any bends in the tubing must be considered in the context of particle losses in the tubing. The selection of the monitoring system should take account of any risk presented by the materials used in manufacturing operation, for examples those involving live organisms or radio pharmaceuticals.***

## **(2) Process Simulation**

A simulation of the process should pinpoint the exact location and height for the probe with the caveat that it does not interfere with the process or get in the operators way. Typically for a BSC or LAF cabinet the sample probe is positioned in the middle of the cabinet since most process zones are in that area and the height of the probe needs to be strategically calculated based on the process. The same applies to critical locations on a filling line where the product is exposed such as when open vials come out of the sterilization oven which are then transported to the filling head and right up to stoppering and capping locations. Probe height is also a critical factor and should be carefully considered in order to get meaningful data from these locations. Do not have the probe too low or too high otherwise your system has too many alarms or none at all. Remember the particle counter function is not to monitor the process but the clean air over the process and to be within a foot of the critical zone.

## **(3) Alarm setting Risk Assessment**

Setting appropriate alarm levels is a process that involves understanding of the activities within the sterile zone. The best approach to understand the risks is to validate the system. This normally occurs during a Performance Qualification (PQ) after the Monitoring System supplier has successfully completed the IQ/OQ. The end user should conduct a proper assessment of the system to understand the limitations and to then set appropriate alarm limits during the PQ. In the case of a filling line or BSC /LAF cabinet the best approach is to setup the process and have settle plate and air sampler also setup as part of the validation (to monitor the process).

Trial runs are conducted under normal operational conditions then separate fail runs are conducted. These fail runs are where the operator purposely generates particles over the process where they can be picked up by the particle counter. The resulting spike in counts are examined as well as the recovery time which in general in a grade A environment with a 0.30-0.45m/s down flow is typically within 1-2 minutes where 0 particles are expected to be recorded by the particle counter. The settle plate and air sampler plates are incubated to back up the expectation. Meaningful and useful data is obtained during this process.

A normal correlation is that very low CFU or even 0 CFU should be expected even with high particle spikes. As long as the integrity of the critical zone remains intact and operators are gowned up correctly. The next phase is to repeat this PQ but with the operator not correctly gowned up with wrists exposed or holes in gloves then repeat the test and results typically yield a resultant correlated CFU return with the particle counter spikes.

The purpose of this test is to understand that there is a much lower probability of CFU counts when operators adhere to proper sterile gown up procedures and when there is a spike in counts. The overall probability of product contamination is nearly negligible when the compounding or sterile activity is also taken into consideration with respect to vial or connection opening and the risk that viable particles enter the vial opening and that the cleanroom is performing under the specified certified conditions. With this PQ validation process end users can establish reasonable and justified alarm limits as they have validated their process and understand the risks and limitations of their facility, personnel, equipment process tools and the process.

#### (4) Setting APPROPRIATE Alarm Limits

The bottom line is that a good Monitoring System should be able to capture transient events as well as repeatable events, interventions and adverse trends. QA Management should be looking for trends rather than one off singular events. The system should immediately notify operators of an action limit event so that they can figure out what is going on at the moment in time as to what caused the particle excursion. Locations integrating video cameras can play-back to the time of particle events to verify the probable root cause, such as user intervention, broken vials or if a trend is noticed then possible failure in the environment containment integrity and HEPA filtration. In the end it is all about root cause investigations and effective Corrective Action Preventative Action (CAPA) using meaningful data. If something goes wrong, you capture that event, you notify the right people immediately and importantly you have enough information to help you investigate the root cause efficiently. This is obvious but it doesn't hurt to restate this because a lot of end users get hung up on setting limits to correspond to the GMP Annex 1 or ISO 14644-1 tables, when in reality they should be set to detect events and adverse trends. See clause 20 below from Annex 1, it DOES NOT SAY set "action and alert limits" per the classification table.

***20. Appropriate alert and action limits should be set for the results of particulate and microbiological monitoring. If these limits are exceeded operating procedures should prescribe corrective actions.***



Fig 7 End User validation (PQ) of remote Particle Counter location in a Cleanroom – Grade B Environment.

It is the combination of the alarming features in the Monitoring System software and SOPs that help to achieve compliance. Refer again to Annex 1 clause 20, the regulators are just as interested in the SOPs that come to life and are executed when the alert or action limits are exceeded. Meaning they are more interested in what you do when an alarm is triggered and they want to see how the reaction is implemented to investigate the issue. ***“If these limits are exceeded operating procedures should prescribe corrective actions”.***

Another issue with following GMP and ISO 14644-1 tables if you report counts per cubic feet then you are expected to have sampled a 1 cubic foot of air, this is no problem with 1 cfm remote particle counters as it means they update a count sample report every 60 seconds. The problem is if data is expected in counts per m<sup>3</sup>, sampling a cubic meter of air takes nearly 36 minutes using a 1 cfm flow rate not so good when starting a process and having to wait 36 minutes to find out the environment is outside alarm limits.

Clause 12 and clause 9 below from PICs GMP Annex 1 clearly releases users from having to take cubic meter samples when monitoring.

***9. For Grade A zones, particle monitoring should be undertaken for the full duration of critical processing, including equipment assembly, except where justified by contaminants in the process that would damage the particle counter or present a hazard, e.g. live organisms and radiological hazards. In such cases monitoring during routine equipment set up operation should be undertaken prior to exposure to the risk. Monitoring during simulated operations should also be performed. The Grade A zone should be monitored at such a frequency and with suitable sample size that all interventions, transient events and any system deterioration would be captured and alarms triggered if alert limits are exceeded. It is accepted that it may not always be possible to demonstrate low levels of >5.0 µm at the point of fill when.....***

***12. The sample sizes taken for monitoring purposes using automated systems will usually be a function of the sampling rate of the system used. It is not necessary for the sample volume to be the same as that used for formal classification of clean rooms and clean air devices.***

There are problems though when using 1 minute samples and reporting counts/ft<sup>3</sup>, especially around 5 µm particle counts if trying to adhere to the Annex 1 or ISO 14644-1 Classification tables. 5µm counts are 20 or 29 respectively for action alarms. If you try to report the data as counts per m<sup>3</sup>, you will always get an action alarm when there is a single 5µm count (1 count at 5µm in a minute sample will give you over 35 counts/m<sup>3</sup> when normalized which is >20 counts). So reporting particle count data for a sixty second sample interval as count/m<sup>3</sup> is not practical or easy to alarm on. Rolling averages are used but they do not address the initial problem when a process is first started. The alarm limits really do need to be set based on a risk assessment and the process baseline results during the PQ.

The following alarm strategy has been implemented to address the above problem. Importantly regulatory authorities will expect the end user to support any alarming strategies using sound scientific evidence, in other words historical particle count data, qualification data and to also implement SOP's of actions once an alarm has been triggered. The example below have been used in numerous Lighthouse Monitoring Systems across the World over several years in manufacturing facilities that go through regulatory audits.

## Setting Limits based on good science and PQ data and using trigger SPC's

This strategy allows for occasional events of a single 5µm particle and was developed based on the historical particle count trend observed over time and PQ observation's using statistical process controls (SPC's). The Strategy shows a good understanding that we are monitoring not classifying. Allowing for the odd 5µm particle count and transient events, it will also capture a 5µm excursion and trend. The 5µm alert limit will notify end user of an adverse trend if one exists before it becomes critical. The 0.5µm alert limits will capture adverse trends based on setting appropriate alarm levels (after PQ). Note the below set-points are starting points and should be reviewed and fine-tuned during the PQ process and also reviewed regularly. If your system never alarms then you may have set the limits too high therefore effectively making the Monitoring System redundant. Remember PICs GMP Clause 20.

***"Appropriate alert and action limits should be set for the results of particulate and microbiological monitoring".***

1. Set alert and action limits on a 60 second sample interval. Expected baseline is zero.
2. Set 5µm action limit > 2 Counts /ft<sup>3</sup>.
3. Set 5µm alert limit >1 Count/ft<sup>3</sup>  
*This is not in agreement with Annex 1 classification table if you were to normalize the data to m<sup>3</sup>. This does however allow for the occasional 5µm particle count per the final sentence in clause 9 of Annex 1.*
4. Set 0.5µm action limit > 100 Counts/ft<sup>3</sup>.
5. Set 0.5µm alert >50 Counts/ft<sup>3</sup> or once historical data has been reviewed.  
*Around the 95th percentile or 2 x standard deviation historical mean counts.*
6. Add statistical process control triggers for X out of Y events, for example this can be 2 out of 3 or 5 out of 7 events, the point is that these SPCs must have some validation around them so you reach a point where X minutes of exceeding alarm limits has a known impact on product quality.  
*Something to cover in the PQ by simulation of real events during processing and the probability of product contamination based on high particle counts versus observed microbial sampling and CFU's at the same location at the same time.*

## Summary

Setting alarm limits in sterile manufacturing processes should not be based from a table which has been established for cleanroom certification/classification. The end user should have a well-defined PQ following a risk assessment approach. Which considers the location of the sample probe, length of sample tubing (to reduce particle loss issues) and monitoring system sample update rate based on the actual particle counter flow rate. All assist to establish a logical baseline from simulated or actual process conditions and have a low probability of risk. In the PQ stage the use of settle plates and active microbial air sampling at the same location as the particle counter location can help yield significant data on the impact of high counts and long trends and the probability of product contamination.

Investing in the time and expertise to address these issues and the development of a well-rounded PQ can literally save thousands of dollars and hundreds of wasted man hours looking for root causes that do not exist in the cleanroom.

Leverage Environmental Monitoring Systems supplier's experiences and their level of knowledge in particle monitoring system installations and setting appropriate alarm limits, sample probe placement with all the regulatory requirements, GMP guidelines and standards that are expected. Selecting a Monitoring System partner should not be an exercise in costs and upfront savings but the true value is in the level of experience they have and the ability of their hardware and software to offer a robust and well supported system as well as their ongoing technical expertise and support.

### Useful References:

- FDA – Sterile Drug Products produced by Aseptic Processing – Current Good Manufacturing Practice (2004).
- FDA – 21CFR11 Electronic Records; Electronic Signatures — Scope and Application (2003)
- PICs GMP Annex 1 Manufacture of Sterile Medicinal Products (2017)
- PICs GMP Annex 11 Computerized Systems (2011)
- PI 011-3 Good Practices for Computerized Systems in GXP” Environments [Pharmaceutical Inspection Cooperation Scheme (PIC/s), 2007]
- GAMP 5 (2008)
- ISO 14644-1 Cleanrooms and associated controlled environments – Part 1 Classification of air cleanliness (1999)
- ISO 14644-1 Cleanrooms and associated controlled environments – Part 1 Classification of air cleanliness (2015)
- ISO 14644-4 Cleanrooms and associated controlled environments – Part 4 Design, construction and start-up (2001)
- PICs – Pharmaceutical Inspection Convention Pharmaceutical Inspection Cooperation Scheme – Guide to Good Manufacturing Practice for Medicinal Products PE-009-13 ( 2017)
- PICs – Pharmaceutical Inspection Convention Pharmaceutical Inspection Cooperation Scheme – PICs Guide to Good Manufacturing Practice for Medicinal Products in Healthcare Establishments PE-010-4 (2014)
- CFR 21 Parts 210 and 211
- The Good Automated Manufacturing Practice (GAMP) Guide for Validation of Automated Systems, GAMP 4 (2001)
- ISO 14971:2002 Medical Devices- Application of risk management to medical devices (2001)

## Biography – Jason Kelly Director of Systems – Lighthouse Worldwide Solutions

20 Years Management positions in Environmental Monitoring Systems Service, Design, Installation, Validation and ongoing support. Has worked on many Projects for top Life-Science companies assisting in procurement, delivery and compliance to ensure regulatory acceptance. Worked across the World on many projects in the UK, Ireland, Europe, Australia and now resides in Oregon USA. He can be contacted by email on [jasonk@golighthouse.com](mailto:jasonk@golighthouse.com) or on LinkedIn and always welcomes queries and questions on Monitoring Systems connected to particle counters or environmental sensors and regulatory compliance.

