





Operators Manual



Lighthouse Worldwide Solutions

ApexZ3 and ApexZ50

Airborne Particle Counters





Operators Manual

Copyright ©2018 by Lighthouse Worldwide Solutions. All rights reserved. No part of this document may be reproduced by any means except as permitted in writing by Lighthouse Worldwide Solutions.

The information contained herein constitutes valuable trade secrets of Lighthouse Worldwide Solutions. You are not permitted to disclose or allow to be disclosed such information except as permitted in writing by Lighthouse Worldwide Solutions.

The information contained herein is subject to change without notice. Lighthouse Worldwide Solutions is not responsible for any damages arising out of your use of the LMS program.

LMS[™], LMS XChange[™], ApexZ[™], ApexZ3[™], ApexZ50[™] are trademarks of Lighthouse Worldwide Solutions. Microsoft[®], Microsoft Windows[™] are trademarks of Microsoft Corporation.

Manufactured by: Lighthouse Worldwide Solutions 300 West Antelope Road White City, Oregon 97503

LWS Part Number: 248083466-1 R11

Table of Contents

Safety	1
About this Manual	1
General Safety	1
Installing the battery	4
Compliance	9
FCC Notice	9
CE Notice	
Introduction	
Overview	
Specifications	
Enclosure and Ports – ApexZ3	
Enclosure and Ports - ApexZ50	
Optional Printer	
LED Handle*	
Operating ApexZ	
Sample Settings Dashboard	
Manual Sample	
Sample Settings	
Channel Settings	
Save Manual Sample	
Start Manual Sample	
Preset Samples	
Assign Preset to Location	
Start Preset Sampling	
Sample Plans	
Select a Sample Plan	
Select a Location	
Save Sample Plan	
Certification Sample	
Quick Certification	
View Options	
Data Table View	

Sample Grid View	. 32
Histogram View	. 35
Graph View	36
Environmental View	37
Sampling with Data Tags	38
Data Records Screen	41
View Data Records	41
Status Bar	43
Data Tags	. 44
Export Data Records	. 45
Make Report from Data Records	. 46
Filter Data Records	. 49
Delete all data records	. 55
Reports	. 56
Report Screen	. 56
Filter Report List	. 58
Export Report	. 60
Delete a Report	. 61
Audit Trail	. 62
Settings	. 65
Instrument Settings	. 67
Network Time Protocol	. 68
Set Current Date	. 68
Set Current Time	. 71
Audio/Visual	. 72
Volume	. 72
Brightness	. 72
Language	. 72
Decimal Mark	72
Digit Grouping	72
Options	73
' Lock Global Channels	
Tool Tips	74
User Options	75

Sample Options	
Zero Count Options	
Data Tag Name Options	
Audit Trail Comments	
Output	
USB	
Print – On-Board Printer	
Print – Network Printer	
CIFS	
Output Formatting - Print	
Output Formattingcsv/.xlsx	
Output Formattingpdf	
Configuration	
Save ApexZ configuration file to USB key	
Load ApexZ configuration file from a USB key	
Connectivity	
LMS XChange Software 2.x	
Network - Ethernet	
Network - Wi-Fi	
Active Directory AD/LDAP	
CIFS -Common Internet File System	
Network Printer	
3 rd Party API - Application Programing Interface	
Locations and Groups	114
View Locations	
Create New Locations	115
Edit a Location	116
Delete a Location	
Add a New Group	
Assign Locations to a new Group	
Edit a Group	
Delete a Group	
Presets	
Add a Preset	

Enter Preset Sample Settings	
Edit Channel	
Enter Warning and Alarm settings	
Environmental Alerts	
Assign Locations to the New Preset	
Preset Name	
Locations linked to Presets	
Sample Plans	
Add New Sample Plan	
Preset - Sample Settings	
Delete Sample Plan	
Certifications	
Add a Certification	
Certification Types	
Edit a Certification	
Delete a Certification	
Users & Permission Levels	
Users	
Add a New User	
Edit User	
Delete a User	
User Login	
User Logout	
Administrator Permission Level	
Add new Permission Level	
Edit Permission Levels	
Assign AD Group to Permission Level	
Prioritize Permission Levels	
Service	
Maintenance	
About Screen	
Service Login	
Calibration	
Trouble Shooting	

Clean and Disinfect	182
Terminology	184
Appendix A: Limited Warranty	186

*Patent Pending

Safety

About this Manual

This manual describes the setup, operation and use of the Lighthouse Worldwide Solutions ApexZ family of Portable Airborne Particle Counters. The word "ApexZ", "unit", or "instrument" may be used in place of a specific model in this ApexZ Operators Manual.

General Safety

Warnings and cautions are used throughout this manual and the reader should become familiar with the meaning of a warning before operating the particle counter. Most warnings will appear next to the subject or step to which it applies. Take care when performing any procedures preceded by or containing a warning. The classifications of warning are defined as follows:

- Laser pertaining to exposure to visible or invisible LASER radiation.
- Electrostatic pertaining to electrostatic discharge.
- Network Connection pertaining to communication ports and instrument damage.

WARNING: There are no user serviceable components inside the particle counter.

Laser Safety

This product is considered to be a Class 1 LASER product (as defined by FDA 21 CFR, §1040.10 and IEC 60825-1:2014) when used under normal operation and maintenance. Performing service on the internal sensor can, however, result in exposure to invisible radiation. The particle counter has been evaluated and tested in accordance with EN 61010-1:2012, "Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use" and IEC 60825-1:2014, "Safety of LASER Products". For further technical assistance, contact our Technical Support Team at 1-800-945-5905 (USA Toll Free), or 1-541-770-5905.



Laser Warning Label



WARNING: There are no user-serviceable components inside the particle counter. The use of controls, adjustments or procedures other than those specified within this manual may result in personal injury and/or damage to this instrument.



Lighthouse Worldwide Solutions

Battery Safety

Battery Type:

National Power Sm-Energy[®] Model: SM221 Lithium-Ion Rechargeable Battery, 14.4V, 6.45Ah, 93Wh



Battery Level Indicator Button

Press here to display blue LED battery charge level 25-50-75-100%



Note: The ApexZ Smart Batteries have a power-saving function where they will go into a hibernation mode after not being used for a period of a few weeks. While in this mode pressing the Test button will not show any charge present on the battery, and the battery will not power the unit. To "wake" the battery, power must be applied to it, such as being inserted into an ApexZ unit that has DC power applied.

ApexZ Operators Manual

Installing the battery

Insert battery 1 into the lower battery slot.

Note: ApexZ has dual battery slots and can operate with either one or two batteries installed. With two batteries installed both batteries will self-adjust and drain their charge evenly together.

Insert battery 2 into the upper battery bay.

You will feel a slight pressure and clicking as the battery seats down fully in the battery bay.





Close the battery door. Tighten down the battery door screw.



Hot swapping the batteries

If there are two batteries installed in the unit, the ApexZ will allow the user to remove one battery while the instrument is powered on and sampling. This allows the user to continue to sample with one battery while recharging the second battery. A charged second battery may then be reinstalled into the instrument while ApexZ is still sampling.

Hot swapping allows the user to install a newly charged battery and then remove a low charge battery while the instrument is sampling under battery power.

Use the optional battery charger (PN 017789753: Dual Battery Charger) to charge batteries outside of the ApexZ instrument.

Batteries inserted into ApexZ will always charge while ApexZ is plugged into the AC power with the AC power adapter and power cord.

WARNING: DO NOT disconnect the power cord or remove the battery while the instrument is powering down. These actions will cause corruption of firmware.





WARNING: Only use batteries provided with the ApexZ by Lighthouse Worldwide Solutions. Use of unauthorized batteries may cause damage to equipment, void warranty, and pose a threat to safety.

Charging the batteries outside of ApexZ

Optional Dual Bay Desktop Charger for charging up to 2 ApexZ batteries at the same time. Batteries not included.




Power Safety

Start Up

Press the power button to turn on the ApexZ particle counter.





Shut Down

To shut down press the power button once. DO NOT HOLD the power button down.



Users are not required to logout before ApexZ is shut down. If Users are enabled, shutting down the ApexZ will logout the current user and will require a user to login when the ApexZ is next started. It is permissible to shut down the ApexZ particle counter from any one of the user interface screens.

WARNING: DO NOT disconnect the power cord or remove the battery while the instrument is powering down. These actions will cause corruption of firmware.

WARNING: Do not press and hold the Power button to shut down. This will cause a hard reset and is not recommended.

WARNING: If the ApexZ is shut down during a sample, that sample will not be completed or saved.

Sampling Safety

WARNING: While the flow path components will not be adversely affected, sampling reactive gasses such as hydrogen and oxygen can be dangerous. Appropriate precautions should be taken. Sampling of any pressurized gasses without the use of properly designed diffuser can cause damage to the instrument and void the warranty. DO NOT allow water, solvents, or other liquids to enter the instrument, as they can damage the instrument and void the warranty.

Operating Safety

WARNING: The use of controls, adjustments or procedures other than those specified within this manual may result in personal injury and/or damage to this instrument. Attempts by untrained personnel to disassemble, alter, modify or adjust the electronics or optics may result in personal injury and damage to the instrument and will void its warranty. There are no user-serviceable components inside the particle counter. Only factory authorized service personnel should repair or service this instrument and its optical system. If replacement of the power supply or its AC power cord is required, replace it only with a power supply or cord having as good as or better ratings than specifications provided by Lighthouse Worldwide Solutions. Attempting to use an underrated power supply or cord can expose the instrument and adjacent equipment to the user to dangerous shock and fire hazards. Failure to meet this requirement will void the CE certification, void the instrument warranty and can result in serious personal injury.

Electrostatic Safety Information

WARNING: Electrostatic discharge (ESD) can damage or destroy electronic components. Therefore, any service or maintenance work should be done at a static-free work station. A static-free work station requires an ESD consultant to evaluate the work environment and propose the equipment and apparel needed for a work station to be successful. Blank Page

Compliance

FCC Notice



This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Class A Digital Devices

NOTE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his expense.

NOTE: Any changes or modifications not expressly approved by Lighthouse Worldwide Solutions could void your authority to operate this equipment.

FCC RF Exposure Statement

This equipment complies with FCC RF exposure limits set forth for an uncontrolled environment. This equipment was verified for RF exposure and found to comply with FCC OET-65 RF exposure requirements.

WARNING: Making changes to the antenna or the equipment is not permitted. Doing so may result in the equipment exceeding the FCC RF exposure guidelines. This equipment must not be co-located or operated in conjunction with any other antenna or radio transmitter.

CE Notice

			C
	EU CE	DECLARATION	I OF CONFORMITY
Manufacture	's Name:		Lighthouse Worldwide Solutions, Inc.
Manufacturer	r's Address:		Lighthouse Worldwide Solutions, Inc. 1221 Disk Drive Medford, OR 97501 USA
Declares the p	product:		
	Product Na Model Nur		ApexZ ApexZ3, ApexZ50
Conforms to t	the Following	; Harmonized Standards:	
	<u>SAFETY:</u>	IEC 61010-1:2010 3 rd Ed.	Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use Part 1: General Requirements
	EMC:	2014/53/EU	Radio Equipment Directive (RED)
			EN 301 489-1, Common technical requirements
			EN 301 489-17, Specific conditions for Broadband Data Transmission Systems
			EN 300 328, Data transmission equipment operating in the 2.4 GHz ISM band and using wide band modulation techniques
		EN 61326-1:2013	Electrical Equipment for Measurement, Control and Laboratory Use
Supplementa	ry informatio	n	
This declaration	on of conforn	nity is issued under the sole	responsibility of the manufacturer.
			0
2017-09-22, N	Medford OR, I	JSA	Jeny Syrach
Date / Locatio	n Issued		Jerry Szpak VP of Engineering

Introduction

Overview

This operating manual introduces the Lighthouse **ApexZ** portable particle counters ApexZ3 and ApexZ50 and includes instructions for inspecting, installing, using and maintaining the instrument.



Description

The ApexZ portable particle counter comes standard with 6 particle-size channels (ApexZ3: 0.3, 0.5, 1.0, 3.0, 5.0, 10.0 microns), (ApexZ50: 0.5, 0.7, 1.0, 3.0, 5.0, 10.0 microns) and a flow rate of (ApexZ3: 1.0 CFM), (ApexZ50: 100LPM). The instrument uses a LASER diode light source and LASER beam shaping optics to illuminate a cross section of the air flow path. As particles move along this path, they cross the LASER beam and scatter light. The scattered light rays from particles cutting through the beam are steered to a photodiode using a system of collection optics. The photodiode converts the image into a current which is converted to a voltage and amplified by an electronic circuit. The amplitude of the voltage pulse is proportional to the light scattered, which in turn is proportional to the size of the particle. The voltage pulses created by the particle are then processed by additional electronics that analyze the height of each pulse and therefore, the size of each corresponding particle. The result is that the number of particles of various sizes is determined. The ApexZ is effective in both ultra-clean areas (such as ISO Class 2 or Grade A) and in more traditional clean zones rated as ISO Class 3 or higher. Refer to Specifications in this manual for additional instrument information.


Accessories Included

ApexZ Includes: USB flash drive with lanyard (includes ApexZ3/ApexZ5 Read Me First and Operators Manual), ISO probe (ApexZ3: 1 CFM, ApexZ50 100 LPM), Zero count filter (ApexZ3 1 CFM fittings, ApexZ50 100 LPM fittings), USB to Micro-B cable, Exhaust Barb, Printed Instrument Read-Me-First, 10ft Bev-A-Line sample tubing, 24 VDC 5.0 amp power adapter and power cord, 1 Li-Ion Smart battery (2 battery slot capacity), 2x rolls of printer paper (units with printers only), Calibration certificate, Software Download Instructions.


Specifications

Size Range	ApexZ3: 0.3 - 10.0 μm.
0.20 1.080	ApexZ50: $0.5 - 10.0 \mu$ m.
Standard 6 Channels	ApexZ3: 0.3, 0.5, 1.0, 3.0, 5.0, 10.0 μm.
	ApexZ50: 0.5, 0.7, 1.0, 3.0, 5.0, 10.0 μm.
Flow rate	ApexZ3: 1.0 CFM (28.3 LPM).
	ApexZ50: 100 LPM.
Laser Source	Laser Diode.
Zero Count Level	< 1 count / 5 minutes (meets specifications of JIS B 9921 and reporting requirements of ISO
	21501-4 Annex C).
Concentration Limits	ApexZ3: 1,000,000 particles / ft3 @ 10% coincidence loss.
	ApexZ50: 300,000 particles / ft3 @ 10% coincidence loss.
Counting Efficiency	ApexZ3: 50% @ 0.3 μm, 100% for particles > 1.5x - 2x first channel particle size.
	ApexZ50 50% @ 0.5 μ m, 100% for particles > 1.5x – 2x first channel particle size.
Calibration	Meets/Exceeds ISO 21501-4 compliance.
Self-Diagnostics	Laser power/current/supply, flow out-of-range, no flow, PA background/supply/health.
Count Modes	Sample, Concentration, Beep.
Data Storage	10,000 sample records, 2500 locations, 250 Groups, 500 Presets, 500 Sample Plans, 500
2444 0101480	Certifications, 10 Permission levels, 100 reports, and 1000 users.
Communication	Ethernet, Wi-Fi (802.11 b/g/n), USB key.
Supporting Software	LMS XChange 2.x
Environmental Sensors	Temperature/Relative Humidity probe: $32-122 \degree F$ (0-50 °C) $\pm 1 \degree F$ (0.5 °C), $15\%-90\% \pm 5\%$ RH.
Touch Screen Display	7.0 inch (8.89 cm), 1200 x 800 IPS color touch screen.
Speakers	2x adjustable speakers up to 70 dB to indicate status and alerts.
Carry Handle	With RGB LED's for status and alarms.
Print	Optional embedded panel-mount printer; IPP (Internet Printing Protocol)
Reports	FED-209E, ISO-14644-1:1999, ISO-14644-1:2015 and EU GMP regulatory pass/fail reports
Reports	displayed on screen or printed.
Alert	User configurable warnings and alarms.
Sample Inlet	ApexZ3: 3/8 inch Inlet nozzle.
campic met	ApexZ50: 1/2 inch inlet nozzle.
Sample Output	Internally filtered to HEPA standard (>99.97% @ 0.3 μ m); with rotatable exhaust fitting for
campic catput	directed flow.
Vacuum Source	Internal blower, automatic flow control.
Enclosure	Chemically Resistant Polycarbonate blend, conductive lining, VHP compatible.
Power	External power supply: 24 VDC, 5 A max draw.
Battery	2x Li-lon battery capacity, removable and rechargeable.
Battery Run Time	ApexZ3: 2 batteries 11 hours nominal use, 50% duty cycle no printing.
	ApexZ50: 2 batteries 8 hours nominal use, 50% duty cycle no printing.
Dimensions	12.76" (I) x 6.23" (w) x 8.88" (h), [32.41 x 15.82 x 22.55 cm].
Temperature Range	(50°F - 104°F) (10°C - 40°C), 20% - 95% relative humidity, non-condensing.
Weight	7.05 lbs. (3.19 kg) without batteries, battery weight 1.35 lbs. (0.61 kg).
0 -	
Languages	English, Spanish, Portuguese, German, French, Italian, Russian, Chinese, Korean, Japanese,

Enclosure and Ports – ApexZ3



Enclosure and Ports - ApexZ50



Optional Printer


Load Printer Paper

Open Printer Door

Pinch the left and right printer door tabs and gently pull open the printer door.

Load Printer Paper

Insert the printer paper roll with the paper feeding from the back of the roll.



Close Printer Door

Close the printer door so that the left and right printer door tabs snap shut.

Warning: Changing the instrument language will update the instrument firmware and the printer. It will take 1-2 minutes for the update to complete. During this time make sure the ApexZ instrument has at least 50% battery life remaining or is connected to the AC adapter and do not turn off the instrument while updating. Failure to follow these instructions could cause the printer to become damaged and unusable.

Paper Feed

The printer paper feed button is available on the ApexZ:

Home Screen, Data Screen and Report Screen.



■ 06/18/2019 2:24:27 pm	Admin 1	Ψ 🔉	₽	i i	(Group 1
🜻 μm			Σ		•	Location 0001
0.3			0		≕, ⊑ ,	Sample Plans
0.5			0		—	
1.0			0		<u>"</u>	Manual
3.0			0		Mode: Cycles: Delay:	Sample 0 / 1 00:00:00
5.0			0		Sample: Hold:	00:01:00 00:00:00
10.0			0		Recs:	9 / 10000
						IDLE
♠ =	a	Fee				•


ApexZ Home Screen


Advance the printer paper

Press the Paper Feed button to advance the paper from the ApexZ embedded printer.

Tear off the printed report

Tear the printed report along the serrated edge of the printer paper opening.





LED Handle*

The ApexZ RGB LED handle lights will light up during and after a sample to indicate warning and alarm status. If the handle LED lights are blinking then the ApexZ is currently sampling. If the handle LED light is a solid color the ApexZ has stopped sampling. The sampling LED color that remains lit after a sample is completed indicates the last recorded warning or alarm condition.


Blinking LED colors indicate ApexZ is actively sampling. (Blue also means alarms are disabled)

- Blinking Blue actively sampling
- Blinking Green actively sampling
- Blinking Orange actively sampling
- Blinking Red actively sampling

no instrument faults, (warnings/alarms are disabled). no instrument faults, no warnings, no alarms. no instrument faults, warning limit has been exceeded. alarm limit exceeded and/or faults detected.

Solid LED color indicates ApexZ is not currently sampling.

- Solid Green sample completed
- Solid Orange sample completed
- Solid Red sample completed










*Patent Pending



Blank Page

Operating ApexZ

This is the Home Screen on the ApexZ airborne particle counter. Create manual or preset sample plans and certifications for up to 2500 unique locations.



Home Screen

Sample Settings Dashboard



Select the manual sample button.

If no group is selected (no group name is highlighted blue) then all locations will display in the location column. The first location will be selected and highlighted blue.

Groups and Locations

Edit the Sample Settings

Press the Manual Sample button to edit the sample mode, cycles, delay, sample time and hold time.

×	Sa	mple Settings	Dashboa	ird	•
Group Name		Location		<u>11</u>	123
Group 1		Location 0001	\uparrow	Ŧ	
Group 2		Location 0002			
Group 3		Location 0003			
		Location 0004		Mode:	Sample
		Location 0005		Cycles: Delay:	1 00:00:00
		Location 0006		Sample: Hold:	00:01:00 00:00:00
		Location 0007			00.00.00
		Location 0008			
		Location 0009	\checkmark		
	3 / 250		250 / 2500		

Sample Settings Dashboard

Manual Sample

Manual Sample Settings

Sample Settings Press each sample setting button to edit:	×	Mar	nual					
Press each sample setting button to eart.	Mode	Sample		Max		e: 0.0028 : 40.7746	8	Sample
Mode	Delay	00:00:00	Enter V	alue			•	Channel / Alerts
Delay	Sample Time	00:01:00	1	2	3	delete		vironmental
Sample Time	Hold	00:00:00						Alerts
Hold Time	Cycles	1	4	5	6			
Cycles #/# Units ft ³ /m ³ /liters	Units	m ³	7	8	9			
Volume	Volume	0.028		0	•	Apply		

Press the Channel/Alerts button to display the Channel Settings screen.


Channel Settings

Press any green channel size button to disable that channel and it will highlight grey and strike through indicating it has been disabled. Press again to enable that channel.

Press any Warning or Alarm threshold button to set those channels respective thresholds. Select Enable or Disable. Enter the Warning or Alarm threshold limits respectively. Press Apply.



×	М	anual		
Channel	Warning (ft ³)	Alarm (ft ³)		Sample
0.3	0.0	0.0		Channel Alerts
0.5	10.0	20.0	2	Environmenta / Alerts
1.0	0.0	0.0	\backslash	/ Alens
3.0	0.0	0.0		
5.0	0.0	0.0		
10.0	0.0	0.0		



Alarm Limit


×	2	0	
	Enable	Disable	
Enter Alarm			
1	2	3	delete
4	5	6	
7	8	9	
	0		Apply

Save Manual Sample Press the green Save button to display the Home screen with the manual sample settings. Manual Ma

Manual Sample - Channel/Alerts

Start Manual Sample

Press the green Start button to start manual sampling when at the selected location.



■ 06/18/2019 Admin 1	· 박 🔉 맛 대 ୩	() •	Group 1 Location 0001
μm	Σ		
0.3	0	=, ₽ 2	Sample Plans
0.5	0		
1.0	0	<u>"</u>	Manual
3.0	0	Mode: Cycles: Delay:	Sample 0 / 1 00:00:00
5.0	0	Sample: Hold:	00:01:00 00:00:00
10.0	0	Recs:	9 / 10000
			IDLE
 ★ ★ 	Feed		•

Home Screen – Manual Sample Loaded

Preset Samples

To run preset sample settings from the Home screen press the Group/Location button to display the Sample Settings Dashboard.

Press the Preset button and it will highlight blue indicating that it has been selected.

All saved preset names will be displayed on the preset name list. Select a group name and only those locations assigned to that group will display on the location list. If no groups are selected then all locations will be displayed.

Locations assigned to a preset will appear with the link icon indicated that a preset is currently assigned to that location.


Saved presets will display in the Preset Name column.

Select a location with no linked preset and all presets will be displayed.

× Sample Settings Dashboard Group Name Location 123 Group 1 Location 0001 \uparrow Preset Name Group 2 Location 0002 Preset 1 Group 3 Location 0003 Preset 2 Location 0004 Preset 3 Location 0005 Preset 4 Location 0006 Preset 5 Location 0007 Location 0008 Location 0009 $\mathbf{\gamma}$ 3/250 250 / 2500 5 / 500

Preset Sample Dashboard

Preset Sample Dashboard



Assign Preset to Location

Select a Preset button and it will highlight blue indicating it has been selected.

Save Preset to Location


Press the green Save button to assign the highlighted preset to the highlighted location.



Start Preset Sampling

The selected Preset Sample Settings will be loaded on the Home screen for the selected location.

Press the Start button to start the sampling preset at the selected location.



Home Screen Ready to Sample



Location linked to Preset

× • Sample Settings Dashboard Group Name Location 123 Group 1 6 Location 0001 \uparrow Preset Name Group 2 Location 0002 S Preset 3 Group 3 Location 0003 œ. Location 0004 Location 0005 Location 0006 Location 0007 Location 0008 Location 0009 $\mathbf{1}$ 3/250 250 / 2500 5/500

Linked Preset


If a location linked to a preset is selected then only that preset name will be displayed on the preset name list. A location with a linked preset cannot be modified from the Sample Settings Dashboard.

Unlink a preset to a location

To unlink a preset from a location, return to the Settings then Location screen and edit the Preset selection.

Sample Plans

To run a Sample Plan from the Home screen press the Sample Plans/Certifications button to display the Sample Plan/Certification Dashboard.



Press the Sample Plan button and it will highlight blue indicating that it has been selected.

The Sample Plan Name list will display all saved sample plans.

When no Sample Plan is selected then all locations will be displayed and selectable.

After a Sample Plan is selected then only those locations assigned to that Sample Plan will be displayed and selectable.

×	S	Sample Plan Dashboard			
	E¥.	Location			
	<u>^</u>	Location 0001	\uparrow		
Sample Plan N	lame	Location 0002			
Sample Plan 1		Location 0003			
Sample Plan 2		Location 0004			
Sample Plan 3		Location 0005		Select a sample plan or certification to view sample	
		Location 0006		settings	
		Location 0007			
		Location 0008			
		Location 0009	\mathbf{V}		
	3 / 500	250	/ 2500		

Sample Plan Dashboard

Select a Sample Plan

Select a Sample Plan and it will highlight blue indicating that Sample Plan has been selected.

Locations assigned to that Sample Plan name will display on the Location list.

Select a Location

Select the desired location to be sampled from the Location list and that location will highlight blue indicating it has been selected.




Press the Save button to load the selected Sample Plan and current Location and return to the Home screen.

•

Sample Plan 1 loaded on the Home Screen

06/19/2019 2:15:33 pm	admin	ቁ 🔉 😲	Of i	() •	Group 1 Location 0001
μm		Σ			
0.3		0		≡,	Sample Plan 1
0.5		0			
1.0		0		<u>"</u>	Manual
3.0		0		Mode: Cycles: Delay:	Sample 0 / 1 00:00:00
5.0		0		Sample: Hold:	00:01:00 00:00:00
10.0		0		Recs:	0 / 10000
					IDLE
♠ 🗉		Feed			•



Selecting a Sample Plan and Location

Certification Sample

To run a Certification Sample from the Home screen press the Sample Plans/Certifications button to display the Sample Plan/Certification Dashboard.

		nome	Screen		
■ 06/19/2019 1:43:02 pm	admin 🌵	× V	04 1	•	Group 1
🗢 µm		Σ		•	Location 0001
0.3		Û		■	Sample Plans
0.5		0			
1.0		0		<u>"</u>	Manual
3.0		0		Mode: Cycles: Delay:	Sample 0 / 1 00:00:00
5.0		0		Sample: Hold:	00:01:00 00:00:00
10.0		0		Recs:	9 / 10000
					IDLE
♠ = 1	e 🔶	Feed			•

Press the Certification button and it will highlight blue indicating it has been selected.

The Certification Name list will display all saved certifications. When no certification is selected all locations will be displayed and selectable.

After a certification is selected then only those locations assigned to that certification name will be displayed and selectable.

× Sample Plan Dashboard				
=, 😜	Location			
	Location 0001	\uparrow		
Certification Name	Location 0002			
Certification 1	Location 0003			
Certification 2	Location 0004			
Certification 3	Location 0005	Select a sample plan or certification to view sample		
	Location 0006	settings		
	Location 0007			
3/5	500 Location 0008			
S	Location 0009	\checkmark		
Quick Certification	2	250 / 2500		

Certification Dashboard

Home Screen


Select a Certification

Select a Certification name and it will highlight blue indicating that it has been selected.

Locations assigned to that certification name will display on the Locations List.

Select a Location

Select the current location to be sampled from the locations list and that location will highlight blue indicating that it has been selected.



Press the green Save button to load the selected certification name and current selected location and return to the Home screen ready to start the certification sample.



Certification Loaded on Home Screen

■ 06/19/2019 2:37:44 pm	admin	ф 🔉	ψ Δ +	Ľ	() •	Group 1 Location 0001
μm 0.3			2)		E*	Certification 1
0.5		()			
1.0		()		C	644-1: 2015 Class 5 Operation
					Cycles: Delay: Sample: Hold: Recs:	0 / 1 00:00:00 00:01:00 00:00:00 0 / 10000
						IDLE
♠ =	e 🔅	Feed				

ApexZ Operators Manual

Sample Plan Dashboard


Quick Certification

Quick Certification is a single input screen for selecting all the required sample parameters for a certification report.


Press the Quick Certification button to display the Quick Cert input screen.

×		Sample Plan Dasht	ooard	1
=	E .	Location		
		Location 0001	\uparrow	
Certification Na	ame	Location 0002		
Certification 1		Location 0003		
Certification 2		Location 0004		
Certification 3		Location 0005		Select a sample plan or certification to view sample
		Location 0006		settings
		Location 0007		
	3 / 500	Location 0008		
	37300	Location 0009	\mathbf{v}	
Quick Certif	ication	250	/ 2500	

Quick Certification Edit Screen



Sample Plan Dashboard



Press the green Save button to Save the Quick Certification to the Home screen.

View Options

While sampling, particle count data may be viewed through multiple view options: Data Table View, Histogram View, Sample Grid View, Histogram View, Graph View and Environmental View. From the Home screen press the View Options button to display the View Options Menu.

			Home Scr	een: Vie	ew Optio	ns Ba	r Displa	ayed
Select View Option		= 06/19/2019 2:43:21 pm ■ JM	admin	₽ ≫	망 마 Σ	ï	() •	Group 1 Location 0001
Sample Grid View Histogram View		0.3			0			Sample Plans
Graph View Environmental View		0.5 1.0			0 0		<u>"</u>	Manual
		₩ 3.0			0		Mode: Cycles: Delay:	Sample 0 / 1 00:00:00
Press a view option bu	tton and the 🚞 🔺	<mark>ار ا</mark>			0		Sample: Hold:	00:01:00 00:00:00
home screen will chang	ge to that view	10.0			0		Recs:	0 / 10000
option and settings.								IDLE
		♠ 🗉 🗄	•	Fee	ł			

Data Table View

	■ 06/19/2019 2:49:55 pm	admin	🜵 🔌 🔱	0t "	() •	Group 1 Location 0001
1	μm		Σ			Location ood i
	0.3		480		=. •	Sample Plans
	0.5		420			
	1.0		360		<u>"</u>	Manual
	3.0		300		Mode: Cycles: Delay:	Sample 1 / 1 00:00:00
	5.0		240		Sample: Hold:	00:01:00 00:00:00
	10.0		180		Recs:	1 / 10000
						IDLE
	↑ ■		Feed			Þ

display the data table view settings.

Press the view options gear icon to

The data table view settings can provide up to three columns of data that can be displayed while sampling.

Data Table View Settings

Enable Column

Select Enable to display that column in the data table view. Select Disable not to display that column.

Display Counts, Warning or Alarm

Select Count to display the particle count, or Select Warning to display particle counts exceeding the warning threshold level, or select Alarm to display particle counts exceeding the alarm threshold level.

Cumulative or Differential

Select Cumulative or Differential to display either cumulative or differential data counts.

Raw or Normalize Units


Select Raw or Norm to display data counts either Raw or Normalized Select units of volume for Normalized particle counts.

Save Data Table View Settings

Press the return icon to save the data table view settings.



The ApexZ Home screen Data View Option will display the data specified in 3 columns.



Data Table View Settings

Home Screen, Data Table View 3 columns enabled

■ 06/19/2019 3:17:15 pm	admin	ቁ 🗴 😗	04 📍	Group 1
🗢 μm	Σ / ft³	Δ / m ³	Σ/L	Location 0001
0.3	480.0	2,118.9	17.0	≕ Sample Plans
0.5	420.0	2,118.9	14.8	
1.0	360.0	2,118.9	12.7	
3.0	300.0	2,118.9	10.6	Mode: Sample Cycles: 1 / 1 Delay: 00:00:00
5.0	240.0	2,118.9	8.5	Sample: 00:01:00 Hold: 00:00:00
10.0	180.0	6,356.8	6.4	Recs: 1 / 10000
				IDLE
♠ ≡ ₿	•	Feed		

Sample Grid View

Select the Sample Grid View options icon from the view options drop down list.

Sample Grid View displays locations assigned to a Sample Plan or a certification with a colored background indicating the sampling status for that location.

There will be no gear icon displayed on the Sample Grid View screen.

06/19/2019 3:31:01 pm ψ U۹ Ū۴ admin () Group 1 • Location 0001 cation 0001 Location 0003 Location 0002 amples: 0 / 1 Samples: 0 / 1 Samples: 0 / 1 =, Sample Plan 1 Manual Mode Sample 0'/ 00:00:00 elav 00:01:00 Samo Hold: 00:00:00 Y Recs: / 10000 IDLE Feed A \triangleright

Sample Grid View Location Button

Each location button will display two lines of text:

Location name	Location 0001
Number of samples taken / total number of sample cycles	Samples: 0 / 1

Sample Grid View Screen

Sample Grid View allows the user to easily select the current sampling location by pressing the matching location name button on the Sample Grid View screen.


When selected, the location button will highlight with a yellow border color // indicating it is the current location to be sampled.

- 3 Loc	6/19/2019 :33:38 pn ation 000 ⁻ nples: 0 /	1 1	admin Locatio Sample			bcation 000	() ()	Locat	Group 1 ion 0001
Sail		'	Sample	5.071		impies. 07	lĺ	Samp	le Plan 1
							Mode: Cycles Delay: Sample Hold: Recs:		Manual Sample 0 / 1 00:00:00 00:01:00 00:00:00 1 / 10000
								IDLE	
n	Ш	<mark>4 </mark>	\$		Feed			J	

Select View Options Sample Grid View

Identify your monitoring point location and select that location on the Sample Grid View screen by pressing that location name button. The selected location button will highlight with a yellow border and the Start button will highlight green indicating the ApexZ instrument is ready to sample at the current location. Make sure your sample tube and sample probe are in place at the selected location point.

Press the green start button to start sampling at the current selected location.



Home Screen – Sample Grid View

After sampling completes then move the ApexZ to the next monitoring point. Make sure your sample tube and probe are in place at the new location and then select the matching current location on the Sample Grid View screen. Press the green start button to start sampling at the new location.

Location Color Status

Yellow border = current location


Green back ground = sample completed with no alarms exceeded.

Blue background = currently sampling.

Note if a sample has multiple cycles then the background color will change after the first cycle is completed if the particle count exceeds the alarm thresholds or the sensor goes into fault.

Grey background = location not sampled </

Home Screen – Sample Grid View



Red background = sample completed and the alarm threshold particle count was exceeded at that location. Note a red background color will also display if the instrument was in fault for its laser, flow or calibration due date.

Sample Grid View Sample Data Display

To view data counts while sampling from Sample Grid View, press the location button that is currently sampling and highlighted blue.

Data counts (from the first two columns of the Data View Settings) will display in a pop up window with white text against a blue background.


≡	06/19/2019 3:48:11 pm	a	admin	ψ	×	С ^р	04	<u>0</u>	()	Group 1
	ocation 0001 amples: 1 / 1		Location Samples				ion 000 les: 0 /		•	Location 0003
			0.3 0.5		Σ / ft 480.0 420.0	3)	۲ 2,	' \ / m³ 118.9 118.9	١ŕ	Sample Plan 1
		3 5	.0 6.0 6.0 0.0		360.0 300.0 240.0 180.0))	2, 2, 2,	118.9 118.9 118.9 356.8	Mode: Cycles: Delay: Sample Hold: Recs:	00:00:00
										SAMPLING
n		<mark>41</mark>	٠):00:5):00:5					•

Home Screen – Sample Grid View

When a Sample Plan is running the user will be prompted to start each sample at each location. The user may select any of the available locations from that Sample Plan to be the next location. When the desired next location has been selected and the sample probe is at the matching location, the user may press Start to begin sampling at that location.

Histogram View

Select Histogram View Option



Select Linear or Logarithmic. Select Cumulative (Σ) or the Differential (Δ). Select Raw or Normalized with units of volume.

If Normalized is selected then select units; cubic feet, cubic meters or liters.



Press the back arrow button to apply the settings and return to the Histogram View.


Note: Histogram View displays real time particle count data as a Histogram.

Histogram View Settings




Graph View

Selecting Graph View Option



Graph View Screen




Graph Settings

Graph View Settings

Select Linear or Logarithmic Select Channel size Select Format Raw or Normalized with units ft³, m³, L.



Press the back button to return to the Graph View screen.



Environmental View

ApexZ can sample and collect temperature and humidity data readings while counting particles during a sample session.

Connect an ApexZ Temperature / Humidity probe and view temperature and humidity data while sampling in real time from the Environmental View Option screen.

From the Home Screen press the View Options button. Then press the Environmental View button.



Press the gear button to display the Environmental View Settings Screen



04/30/2021 11.20:17 am ψ Ur admin N Ē14 • Location 0007 Гуре Value perature 23.1 °C Sample Plans 83 tive 42.0 % idity h. Manual Sample 00:00:00 00:00:10 00:00:00 8 2/10000 IDLE

Select Environmental View Option

Environmental View Screen

Feed

04/30/2021 11:13:37 am	admin	ψ 🔌	ው 🗗 🗸	•	Location 0007
🔫 🍄 Type		Value			
Temperature		23.1 °C		=	Sample Plans
Relative		40.8 %			campie i iane
Humidity				<u></u>	Manual
				Mode: Cycles: Delay: Sample Hold: Recs:	Manual Sample 0 / 1 00:00:00 : 00:00:10 00:00:00 2 / 10000
					IDLE
_ =		Feed			

Environmental View Settings Screen

Environmental View Settings


Press the Enable button and it will highlight blue indicating the sensor type has been enabled.

Press Value to enable selection of °F or °C. Select any of the four Alarm or Warning buttons to enable display of those values.

Relative Humidity will be displayed as %.



Press the Back button to return to the Home page Environmental View.



Sampling with Data Tags

See Instrument–Options–Data Tags to create and enable data tag prompts pre-sample or post-sample.

Pre-Sample Data Tag Input (*created and enabled from the Instrument-Options-Data Tags screen*) will display on the ApexZ screen immediately after pressing START.



Select the data tag name to enter or scan.


Scan or touch to input the data tag.

No Scanner:

If the user selects a Saved Data Tag Annotation (created when setting up the Data Tag) the ApexZ will display the Annotation text in the display area. Alternately, the user can touch the input button to enter the Enter Value input screen and the user can input the data tag value manually and press the green Apply button to apply the data tag value.



Once scanned or manually input, press the Save button then the green Start button to begin the Sample.



Post-Sample Data Tag Input (*created and enabled from the Instrument-Options-Data Tags screen*) will be included with the recorded data. During creation of the Post Sample Tags the user can opt to have a popup display on the screen immediately after sampling is complete. Otherwise, the user will have to go to the data screen and input the value there.

Scan or touch to input the data tag.

No Scanner:

If the user selects a Saved Data Tag Annotation (created when setting up the Data Tag) the ApexZ will display the Annotation text in the display area. Alternately, the user can touch the input button to enter the Enter Value input screen and the user can input the data tag value manually and press the green Apply button to apply the data tag value.

×											
Enter Val	ue										
		w	е	r	t	у	u	i		Р	
		а	s	d	f	g h	j	k	I		
	Shift	z	×	c	v	b	n	m		Delete	
	123			Space Apply							



Once scanned or manually input, select the green Save button to save the record and return to the Main screen.

Blank Page

Data Records Screen

View Data Records

The ApexZ can store up to 10,000 sample data records and display them on the data screen. The user can scroll up and down, left and right and view all saved data records. Individual or multiple data records can be selected and a Data Tag can be entered and saved in the data record. The columns displayed in the Data Records screen can be chosen from 30 options, and the display order can be re-organized by the user.

Data Records screen columns displayable

Date and Time stamp (always shown) Sample Name (3 hidden columns) Sample Time Group Location Channels (up to 6 hidden columns displaying particle counts by Channel Size) Environmental (Analog Temp and RH) Alerts (Alarm and Warning conditions) Faults (Flow, Laser, and Service) User Data Tag(s)

Note: Hidden columns can be viewed by clicking on the column expansion arrow buttons to the right of each column heading.

Data Records Screen

3	Sample 🗲				Filter
Timestamp	Name	Sample Time	Group	Loci	Time
04/30/2021, 12:43:41 pm	Manual	00:00:10	n/s	↑)
04/30/2021, 12:43:17 pm	Manual	00:00:10	n/a		Grou
04/30/2021, 12:42:45 pm	Manual	00:00:10	n/a		a 114 17.555
04/30/2021, 12:42:26 pm	Manual	00:00:10	n/a	•	Location
04/30/2021, 12:42:06 pm	Manual	00:00:10	n/a	8	Sample Plar
04/30/2021, 12:41:46 pm	Manual	00:00:10	n/a		
04/30/2021, 12:41:26 pm	Manual	00:00:10	n/s	л. 🏷	Tag
04/20/2021 12:1E-49 nm		00-00-10	`د	¥ .	
ving: 9 / 9	selected: 0		data buffer: 9 / 1	0000	Alert / Faul
♠ 🔳 🔒		eed All	🕒 + 🖻		B

Note: Data records with a strikethrough present across the entire record indicate that the unit's Calibration Due Date had passed when the sample was taken.

Data Record Screen Column Display

The user is able to select which data fields to display on the Records screen and what order the columns are displayed on the screen. Press the gear icon in the upper left corner of the Data screen to open the Show Record Columns screen.

Select the desired columns to display on the left side to add them to the "selected" list. Active columns will be highlighted in blue.

Once columns are added to the "selected" list, the user can adjust the display order by doing a "press-and-hold" on the at the beginning of the line, then dragging the column into the desired placement. Once all columns in the list are ordered properly, press the green Save button to save and return to the Data Records screen.

× Sł	Show Record Columns						
Column : 31		selected : 11					
Column		Visible Columns					
Sample	\uparrow	IIII Sample	1				
Certification		iiii Group					
Sample Plan		IIII Location					
➡ Preset		IIII Channels					
Sample Time		IIII Environmental					
Group		iiii Alerts					
Location		iiii Faults					
Channels	\mathbf{v}	IIII User	Y				

Channels column expanded

Press the right arrow button to expand and the left arrow button to collapse the channel columns.

Slide the horizontal scroll bar or press the right and left arrows to view additional data columns.

Slide the vertical scroll bar or click on the arrows up and down to scroll through data records saved on the ApexZ.

Data screen displaying channel columns expanded

ŝ		Channe -						Filter
	Timestamp	0.3	0.5	1.0	3.0	5.0		Time
	04/30/2021, 12:43:41 pm	80	70	60	50	\uparrow	5	20003
_	04/30/2021, 12:43:17 pm	80	70	60	50			Group
	04/30/2021, 12:42:45 pm	80	70	60	50			
	04/30/2021, 12:42:26 pm	80	70	60	50		•	Locatior
	04/30/2021, 12:42:06 pm	80	70	60	50		0 5	Sample Plar
	04/30/2021, 12:41:46 pm	80	70	60	50		8	dentifik attechtiken
	04/30/2021, 12:41:26 pm	80	70	60	50		0	Taç
	04/20/2021 12-15-48 pm	1	70		50	*		
ew	ing: 9 / 9	selected: 0			data buffer: 9	/ 10000		Alert / Faul
ł		•	Feed All		+ 🗈			G

Environmental, Alerts, and Faults Columns expanded

Press the right arrow button to expand and the left arrow button to collapse the data columns.

Environmental Columns

Alerts Columns

Faults Columns

None None None None None None None

								_
vironmental 🔸		Alerts 🔸				Faults 🔸		
MP	RH	Alarm	Warning	Environmental Alarm	Environmental Warning	Flow	Laser	Serv
22.6 °C	42.0 %	None	None	None	None	None	None	
22.6 °C	42.0 %	None	None	None	None	None	None	
22.6 °C	42.0 %	None	None	None	None	None	None	
23.1 °C	42.0 %	None	None	None	None	None	None	
23.1 °C	42.0 %	None	None	None	None	None	None	
23.1 °C	42.0 %	None	None	None	None	None	None	
23.1 °C	42.0 %	None	None	None	None	None	None	

Use the horizontal scroll bar to view additional data columns to the right or swipe your finger on screen to scroll up or across and display additional rows and columns.

> Note: Swipe the screen or slide the Elevator Bars to view additional columns to the right or left.

Status Bar

Viewing: #/#

Filtered data records / total data records.

Selected:

This represents total data records currently selected / highlighted blue.

Data buffer: #/10,000

This represents total # of saved data records / total capacity (10,000 records).

viewing: 35 / 35

		Sample •			Channels	Filter
	Timestamp	Name	Group	Location	.0	🛁 Time
	06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	\uparrow	
	06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fill 1		Group
/	06/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		
	06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1		• Location
	06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1		Sample Plan
	06/21/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1		Δ .
ds /	06/21/2019, 9:37:55 am	Sample Plan 1	Group 1	Fill 1		🔁 Tag
	06/21/2019, 9:37:45 am	Sample Plan 1	Groun 1	Fill 1	\rightarrow	Alert / Fault
	06/21/2010_0:27:25 am viewing: 35 / 35	select	led: 0	data buffe	: 35 / 10000	
	♠ 💻 - 1	i	Feed All	🔨 + 🗎		Ð
			i	/		
	**		l.			
	selected:	0		🎽 data bu	ffer: 35	/ 10000

Data Screen

Select Data Records

Press anywhere on a data record row and that data record will highlight blue indicating it is currently selected. Select additional rows and they will highlight blue indicating they are also selected now as a group.



To unselect press again and the data record will return to grey/dark grey highlight indicating the data record is no longer selected.

Data records selected

		Sam	ple 🗖			Channels		Filter
	Timestamp		Name	Group	Location	0.3	-	Time
06	/21/2019, 9:39:06 a	am	Sample Plan 1	Group 1	Fill 1	\uparrow	D	
06	/21/2019, 9:38:56 a	am	Sample Plan 1	Group 1	Fill 1			Group
> 06	/21/2019, 9:38:46 a	am	Sample Plan 1	Group 1	Fill 1		-	
> 06	/21/2019, 9:38:36 a	am	Sample Plan 1	Group 1	Fill 1		?	Location
06	/21/2019, 9:38:26 a	am	Sample Plan 1	Group 1	Fill 1		펳	Sample Plan
06	/21/2019, 9:38:05 a	am	Sample Plan 1	Group 1	Fill 1		Δ	
06	/21/2019, 9:37:55 a	am	Sample Plan 1	Group 1	Fill 1	↑	•	Tag
06	/21/2019, 9:37:45 a	am 🖌	Sample Plan 1	Groun 1	Fill 1	\mathbf{v}		
vie	wing: 35 / 35		select	ted: 4	data buffer	→ : 35 / 10000	A	Alert / Fault
	♠ 🔳	<mark>41</mark>	٠	Feed All	💊 + 🖹			₽

viewing: 35 / 35

selected: 4 ----

data buffer: 35 / 10000

Data Tags

Data Tags can be used to add process references to your sample data record. For example data tags "Batch #" and "Lot #" can be created and prompted for either Pre-Sample or Post-Sample.

Data Tags can be entered by a standard USB barcode reader, input via the use of Annotations, or input from the ApexZ touch screen keyboard. Data Tag input is saved with the data record and displayed in columns on the ApexZ Data Records Screen (if the individual Data Tag column(s) are selected to be shown).

Data Tags can also be added to records in the Data Records screen.

Data Records Screen

Select Data Records

Select a single data record or a group of data records and they will highlight blue indicating these records are currently selected.

Press the Tag button to display the Data Tags screen.



	Sample 🔹			Channels I	Filter
Timestamp	Name	Group	Location	0.:	Time
06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	个 【	9
06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fill 1		Group
06/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1	•	Location
06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1		Sample Plan
06/21/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:37:55 am	Sample Plan 1	Group 1	Fill 1	\	Tag
06/21/2019, 9:37:45 am	Sample Plan 1	Group 1	Fill 1	\rightarrow	
06/21/2019_0-37-35 am viewing: 35 / 35	selec	ted: 4	data buffer	: 35 / 10000	Alert / Fault
↑ ■		Feed All	+=		Đ

Add Data Tags to Data Records

Select the data tag to input. -

Scan or touch to input the data tag.

No Scanner:

If the user selects a Saved Data Tag Annotation (created when setting up the Data Tag) the ApexZ will display the Annotation text in the display area. Alternately, the user can touch the input button to enter the Enter Value input screen and the user can input the data tag value manually and press the green Apply button to apply the data tag value.

Once the values are input, press the green Save button to add the Data Tags to the record(s)

×	Data Tags	
📎 test tag		Press here to update
Satch1		Press here to update
🇞 post tag	→ 	Press here to update

×										
Inter Val	ue									
			•	r .	t	у	u	÷.		P
		a		d	1	a h	i	k	1	
	Shift	z	×	c	×	ь	n	-		Delete
	123			Space						iply

Export Data Records

Data records on the ApexZ can be output to the onboard printer (on units with the optional embedded printer), a network-connected printer, a USB key or a CIFS (Common Internet File System). Data exported to a USB and CIFS can be saved as a .csv, .xlsx, .lsdx, or .pdf file, dependent on the settings in the Instrument \rightarrow Output screen.

Select Data Records


Select a single data record or a group of data records.

To export All Data Records just press the export button without selecting any record(s).

	Sample •			Channels	Filter
Timestamp	Name	Group	Location	0.:	ı Tim
06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	^ [™])
06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fill 1		Grou
)6/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1	•	Locatio
06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1	8	Sample Pla
06/21/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1	<u> </u>	
)6/21/2019, 9:37:55 am	Sample Plan 1	Group 1	Fill 1	. 🕚	, Ta
06/21/2019, 9:37:45 am	Sample Plan 1	Groun 1	Fill 1	↓	
eving: 35 / 35	selected	4	data buffer	→ <u>·</u> : 35 / 10000	Alert / Fau
☆ ≡ 1		Feed All	🍾 + 🗎		B

Data Records Screen

Press the Export button to display the Export Screen.



Use this screen to enter a desired file name for the data to be exported then press the Apply

button. Alternatively, press the Skip button to title the file with the default system name.

Name Export File


When exporting data files to a USB or to CIFS, the user has the option to give the file a name. For those output destinations, after pressing the export button the Enter File Name screen will appear.

Enter File Name Screen

Make Report from Data Records


Reports may be created from selecting a continuous set of data records from the Data Screen and then pressing the + Report button. Up to 120 continuous data records may be selected at one time to produce a Data Report. If over 120 data records are selected then the + Report button will become disabled.

After the + Report has been pressed, the user will be asked to set report sample parameters and create a new Report based on the selected data record samples collected. Follow these steps to select the data records, report type, channels, enter a new report name, and save the report.



Note: Selecting data records from a filtered display of data records will not filter the data set included in the generated report. All data records from the oldest to the newest data record set selected will be included in the generated report.


Report Type



Select the Report Type:

ISO 14644-1: 2015 EU GMP ISO 14644-1: 1999 FED 209E (Federal Standard 209E)

Press the Channel button next to display the Report Channels Screen.



Report Channels

ApexZ will enable the required channels and will set the target limits for the chosen report type standard.

Enabled channels will display with a green check box.

Disabled channels will be displayed with the null symbol, a circle with a line through it. These particle channel sizes are not within the certification standard and may not be enabled.

Press the Name button to display the Report Name Screen.

New Report Channels

×			New Report		
	Size	Enable	Target Limit / L	E	Туре
	0.3 µm	\oslash			Channel
	0.5 µm	~	3520	•	Channel
	1.0 µm	\oslash		• -	Name
	3.0 µm	\oslash			T
	5.0 µm	~	20		
	10.0 µm	\oslash			

Report Name

Enter the Report Name from the keyboard. Press the green Apply button to apply report name. Press the green Save button to save the report and return to the Data Screen.



The new report will be saved on the Report screen.



Report Name Screen

Filter Data Records

From the Data Screen:

Press the Time button to display the Date & Time Screen Filter



	Sample			Channels	Filter
Timestamp	Name	Group	Location	0.:	Time
06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	↑ [™]	7
06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fiii 1	1	Group
06/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1	/ ?	Location
06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1	8	Sample Plan
06/21/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1	<u>^</u>	
06/21/201 9, 9:37:55 am	Sample Plan 1	Group 1	Fill 1	↓	Тад
06/21/2019, 9:37:45 am	Sample Plan 1	Group 1	Fill 1	\mathbf{v}	
06/21/2019 0:37:35 om viewing: 35 / 35	Select	ted: 0	data buffer	→ <u>·</u> : 35 / 10000	Alert / Fault
♠ ■	.	Feed All	🐤 + 🗎		Ð


Data Screen

Filter Start date and time

Press the Start Date button and select the starting ______ day, month and year for the filter period using the keypad.

Press the **Start Time** button and enter the start time on the keypad.

Filter Data – Date Keypad



Filter End date and time


Press the End Date button and enter the end date on the keypad.

Press the End Time button and enter the end time on the keypad.

Press the green Apply button and press the green Save button to apply the Date and Time filter to the data records on the Data Screen.







Filter Groups

Press the Group button to display the Filter Group Screen.

	Sample •			Channels I	Filter
Timestamp	Name	Group	Location		Time
06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	个 ¹⁰	
06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fill 1		Group
06/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		1
06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1	•	Location
06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1	8	Sample Plan
06/21/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1	<u> </u>	
06/21/2019, 9:37:55 am	Sample Plan 1	Group 1	Fill 1	`	Тад
06/21/2019, 9:37:45 am	Sample Plan 1	Group 1	Fill 1	\rightarrow	
06/21/2019_9:37:35 om viewing: 35 / 35	selec	ted: 0	data buffer		Alert / Fault
↑ ■		📑	🍾 🕂 🖹		B

Data Screen

Select Group to Filter

Select the Group name(s) to be included in the data records filter and press the adjacent include button. The button will highlight green with a check mark indicating the group name has been included in the filter.


Press the include button a second time to remove selection. A grey back ground color indicates no filter on this parameter.

Press the green Save button to save the filter settings and return to the Data Screen.

Group Filter

×		Group Filter	
	Reset Fill	tter 👌	
	Include	Group Name	
		All locations	
	×	Group 001	
		Group 002	
		Group 003	

Filter Locations



Filter Sample Plans and Certifications

Press the Sample Plan button to display the Filter Sample Plans screen.



Data Screen

Filter Sample Plans and Certifications



The button will highlight green with a check mark indicating it will be included in the data record filter.

Press the green Save button to activate the filter settings and return to the Data Screen.



Reset Filter 👌 Reset Filter	
Include Sample Plan Name Include Certification Name	
Filling Line Certification 1	
Sample Plan 1 Certification 2	
Sample Plan 2 Certification 3	
Sample Plan 3	

Filter Tags

Press the Tag button to display the Tag Filter screen.

Data Screen

	Sample 🔹			Channels	Filter
Timestamp	Name	Group	Location		Time
06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	↑ [™])
06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fiii 1		Group
06/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1	•	Location
06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1	8	Sample Plan
06/24/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1	<u> </u>	
06/21/2019, 9:37:55 am	Sample Plan 1	Group 1	Fill 1		Tag
06/21/2019, 9:37:45 am	Sample Plan 1	Group 1	Fill 1	\rightarrow	
06/21/2010 9:27:25 pm viewing: 35 / 35	selec	ted: 0	data buffer	: 35 / 10000	Alert / Fault
↑ ■		Feed All	💊 + 🗎		G

Tag Filter Screen

Select the Include button next to the Tag name(s) to be included in the data records	× Reset F		Filter	•
filter.	Reser			
	Include	e Data Tag Name	Value	
		Batch	137900888	
		Batch	5678	
The button will highlight green with a		Comment	ок	
check mark indicating it will be included in	×	Comment	Alarm ack	
the data record filter.		Notes		
		Sample ID	2577388	
Press the green Save button to activate	 ✓ 	Sample ID	1234	
the filter settings and return to the Data				
Screen.				
6				


Filter Alert/Fault

Press the Alert/Fault button to display the Alert/Fault Screen.

	Sample 🕒			Channels I	Filter
Timestamp	Name	Group	Location		Time
06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	^ [™]	
06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fiii 1		Group
06/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1	•	Location
06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1	8	Sample Plan
06/21/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:37:55 am	Sample Plan 1	Group 1	Fill 1	L. 🔪	Tag
06/21/2019, 9:37:45 am	Sample Plan 1	Groun 1	Fill 1	\rightarrow	
06/21/2019_9:27:25 am viewing: 35 / 35	selec	ted: 0	data buffer	: 35 / 10000	Alert / Fault
↑ ≡		📑	🍾 + 🖹		B

Data Screen

Filter Alert / Fault








Note: Press the Reset Filter button to disable all included filters.

Delete all data records

ApexZ will store up to 10,000 sample data records and after the database is full then the newest data records will copy over the oldest data records, wrapping the database.


Data records may not be deleted individually.



	Sample •			Channels K	Filter
Timestamp	Name	Group	Location		Time
06/21/2019, 9:39:06 am	Sample Plan 1	Group 1	Fill 1	↑ [™])
06/21/2019, 9:38:56 am	Sample Plan 1	Group 1	Fiii 1		Group
06/21/2019, 9:38:46 am	Sample Plan 1	Group 1	Fill 1		
06/21/2019, 9:38:36 am	Sample Plan 1	Group 1	Fill 1	•	Location
06/21/2019, 9:38:26 am	Sample Plan 1	Group 1	Fill 1	8	Sample Plan
06/21/2019, 9:38:05 am	Sample Plan 1	Group 1	Fill 1	<u>^</u>	
06/21/2019, 9:37:55 am	Sample Plan 1	Group 1	Fill 1	``	Тад
06/21/2019, 9:37:45 am	Sample Plan 1	Group 1	Fill 1	\rightarrow	
06/21/2010 9:27:25 pm viewing: 35 / 35	select	ted: 0	data buffer	: 35 / 10000	Alert / Fault
↑ ■		Feed All	💊 + 🗎		Đ

Data Screen

- - Confirm Data Records Delete




Reports

Report Screen

The Reports screen displays a list of all saved certification reports.

Reports that have Passed their certification conditions will be highlighted green and display "PASS" in the Results column. Reports that have failed their certification conditions will be highlighted red and display "FAIL", indicating they have failed.



Expand Result Columns

Press the arrow button to expand the Result _ columns.

Use the horizontal slider button or arrow buttons to display more columns to the right or left.

Result Columns Expanded



Report data columns displayed include date created, report name, and fault results for over limit data counts, minimum sample volume not met, minimum samples not taken, minimum sample time not sampled, and instrument fault.

Report Result Column Expanded

Over Limit, Minimum Volume, Minimum Locations, Minimum Samples, Min. Sample Time, Instrument Fault


	Result 🖸	$\langle \rangle$				
Created	Overlimit Fault	Min Volume Fault	Min Loc Fault	Min Samples Fault	Min Sample Time Fault	Instrument Fault
06/28/2019, 1:07:55 pm	Fault	None	None	None	None	None
06/28/2019, 1:00:06 pm	Fault	None	None	None	None	None
06/27/2019, 2:24:20 pm	None	None	None	None	None	None

Filter Report List

The list of reports displayed on the Reports screen can be filtered by the date the report was created, the report type or the report pass/fail status.

Filter report date and time

From the Reports Screen press the Created button to go to the Filter Created Screen





Reports Screen

Select date, month and year.

Select the Start/End date button and the date keypad will be available to enter the start/end dates for the desired data records filter.

Keset Filter × Mon Tue Wed Tuu Fri Sat Sun Start 30 01 02 03 04 05 05/07/2019 1:26:38 pm 06 07 08 09 10 11 12 End 13 14 15 16 17 18 19 05/07/2019 1:38:27 pm 20 21 22 23 24 25 26 07 08 09 30 31 01 02 03 04 05	×	× Created Filter								
Start 2> 30 01 02 03 04 05 05/07/2019 1:26:38 pm 06 07 08 09 10 11 12 End 13 14 15 16 17 18 19 05/07/2019 1:38:27 pm 20 21 22 23 24 25 26 27 28 29 30 31 01 02	Reset Filter		<		м	ay 201	9		>	
End 13 14 15 16 17 18 19 05/07/2019 1:38:27 pm 20 21 22 23 24 25 26 27 28 29 30 31 01 02	Start									
End 20 21 22 23 24 25 26 05/07/2019 1:38:27 pm 27 28 29 30 31 01 02	05/07/2019	1:26:38 pm	06	07	08	09	10	11	12	
05/07/2019 1:38:27 pm 27 28 29 30 31 01 02	End		13	14	15	16	17	18	19	
27 28 29 30 31 01 02	05/07/2019	1:38:27 pm	20	21	22	23	24	25	26	
03 04 05 06 07 08 09			27	28	29	30	31			

Select the Start time and End time buttons and use the Time keypad to enter the report filter start and end times. Select am/pm. Press the Apply button.

Press the Save button to return to the Reports screen.



Select Time

-

Filter Certification Type

Click to Enable the certification type(s) to be included in the filter. The button will highlight green with a yellow check mark indicating that certification Type is included in the data records filter.

Press the Save button to apply the filter and return to the Reports screen.



Filter Pass/Fail

Click to Enable the Pass/Fail Type(s) to be included in the filter; Cannot Classify, Fail or Pass.

The button will highlight green with a yellow check mark indicating that type is included in the Report list filter.

Press the button a second time and the button will change to a grey color indicating it is not included in the data records filter.

Press the green Save button to apply the filters and return to the Reports Screen.



If you do not press the Save button then no new filters will be applied. If you press the red x button to cancel then no filters will be applied. If no Types are included then no filters will be applied and all Reports will be displayed on the Reports screen.



Note: Press the Reset Filter button to disable all included filters.

Select Certification Type Filter


×	Type Filter									
	Reset Fil	Reset Filter 👏								
	Include	Туре								
	~	EU GMP								
		FED 209E								
		ISO 14644-1: 1999								
		ISO 14644-1: 2015								

Select Pass/Fail Filter

×		Pass / Fail Filter	
	Reset Fil	ter 👌	
	Include	Туре	
		CANNOT CLASSIFY	
	\checkmark	FAIL	
		PASS	
			i i

Export Report

Select Report to Export



Select the Export Type: Printer, Network Printer, USB, or CIFS. File format and setup for the USB and CIFS options are set in the Settings -> Instrument \rightarrow Output screen.


```
Press the Confirm Export button to start the export.
```

Export Output

Name Export File

When exporting data files to a USB or to CIFS, the user has the option to give the file a name. For those output destinations, after pressing the export button the Enter File Name screen will appear.

Use this screen to enter a desired file name for the data to be exported then press the Apply button. Alternatively, press the Skip button to title the file with the default system name.



Enter File Name Screen



Delete a Report

Select a report to delete and it will highlight blue indicating the Report has been selected.

Press the Delete Report button to display the

Delete report confirmation screen.

Select Report to Delete


Result • Filter Created Pass / Fail Туре Name Ê 06/28/2019, 1:07:55 pm ISO 14644-1: 2015 ISO 14644-1: 2015 06/28/2019, 1:00:06 pm Cert 3 FAIL Create 06/27/2019, 2:24:20 pm Cert 2 PASS EU GMP Туре 06/27/2019, 2:17:49 pm Cert 1 FAIL ISO 14644-1: 2015 Pass / Fai ewing: 4 / 4 data buffer: 4 / 100 Feed G

Reports can either be deleted one at a time using the Delete Report button, or all at once using the Delete All button. If no reports are selected then the delete report button will be disabled.

Press the green check box to delete the named Report. Or press the red x box to cancel with no report deleted.

Press the Delete All button to delete all Reports listed in the Reports screen. Press the green check box to delete all reports.

Confirm Report to Delete





Audit Trail Cor

 \mathbf{v}

 \rightarrow

All: 149 / 20000

Filter

G

Audit Trail Screen

Printing report Room1 successful by admin.

USB .xlsx record export successful by admin

Audit Trail

The Audit Trail Screen creates a timestamp and message for every user action and data record recorded on the ApexZ instrument. Every time a sample is recorded or a user changes a current setting or modifies a configuration parameter; a timestamped message will be recorded in the Audit Trail. After 20,000 records have been recorded in the audit trail the next record will replace the oldest record.

The user has the ability to add comments to audit trail messages, filter the displayed messages based on Time parameters, and output the displayed Audit Trail to a USB or CIFS.

Timestamp

05/03/2021, 7:48:56 am

05/03/2021, 7:05:47 am

wing: 149 / 20000

Message

05/03/2021, 7:14:24 am Data report Room1 created by admin.

Press the Reports Screen Button

From the Home screen press the Reports button and then press the Audit Trail button to display the Audit Trail Screen.

Press the side arrows or slide the scroll bar up and down to view the audit trial history, and use the bottom arrows or slide bar to completely view the Audit Trail Comment column

Adding Comments to Audit Trail messages

Highlight a message in the Audit Trail by pressing it.

Press the Comment button to bring up the Audit Trail Comment box. Press the Enter Comment button to bring up the Enter Audit Trail Comment keyboard, or select a pre-saved comment.

×											*
Enter Au	dit Trail Co	omment							u		
	q	w	θ	r	t	У	u	i		р	
		a	s	d	1	9	h j	k	1		
	Shift	z	×	c	v	ь	n	m		Delete	
	123			Space					Ap	oply	

×	Audit Trail Comment
	Enter Comment
Saved A	udit Trail Comments
	No Data

Enter the text desired for the Comment, then press Apply to add the Comment to the Audit message.

Adding Comments to Audit Trail messages also generates another entry in the Audit Trail, which indicates timestamp of the Commented message, the text of the Comment, and the user who created the Comment.

Export Audit Trail to USB or CIFS.

Press the Export button to export the Audit Trail. All messages currently listed in the Audit Trail are exported each time (See Audit Trail Time filtering on next page).

Select export to USB or CIFS.

The format of the exported file will depend on the file format options as set in the Instrument \rightarrow Output screens.



Press the Export button.

Name Export File

After pressing the export button the Enter File Name screen will appear.

Use this screen to enter a desired file name for the data to be exported then press the Apply button. Alternatively, press the Skip button to title the file with the default system name.





Export Audit Trail




If a USB key is inserted then a pop up message will display saying "Preparing Data" and "USB Export" each will count up from 0 to 100%.

Audit Trail Time Filtering

The user is able to filter the messages displayed in the Audit Trail screen in regards to time.

Press the Time button in the Audit Trail screen.

The Time Filter screen will appear. Here the user is able to select Start and End values for Date and Time.




With the Date button selected on the left side, the user can use the calendar shown to select a date.

ad to enter a on to indicate nputs are no 00:00 to 12:5	the par			-	0										
nputs are no	•		uav u	esired			Reset Filter 👌				A	pril 20	21		>
		d, or	•		es		Start		2		Wed 31	Thu 01	Fri 02	Sat 03	Sun 04
	9.59)						04/09/2021	10.11.18 am	o	5 06	07	08	09	10	11
							End		1	2 13	14	15	16	17	18
Ti	me Filter					•	05/03/2021	7-54-07 om	1	20	21	22	23	24	25
	Time						00/03/2021	7.34.27 am	2	3 27	28	29	30	-01	- 02.
↓		10 102 <u>01</u> 2010	00:0	0:00						04					
10:11:18 am	Enter Sta	rt Time I													
	1		3	delete											
7:54:27 am	4	5	6												
	7	8	9												
		0		Apply											
	10:11:18 am	Enter Sta 10:11:18 am 1 7:54:27 am	Time 10:11:18 am Enter Start Time 1 2 4 5 7:54:27 am 7	Time 10:11:18 am 7:54:27 am	Time 00:0:00 10:11:18 am 1 2 3 delete 1 2 3 delete 4 5 6 7 7:54:27 am 7 8 9 8000000000000000000000000000000000000	Time 00:000 10:11:18 am 1 2 3 delete 1 2 3 delete 4 5 6 7 7:54:27 am 7 8 9 1000000000000000000000000000000000000	Time 10.11:18 am 7:54:27 am	Time Filter 05/03/2021 10:11:18 am 01:00:00:00 1 2 3 delete 4 5 6 0 7:54:27 am 7 8 9 0	Time Filter 05/03/2021 7:54:27 am 10:11:18 am Image: Control of the second	Time Filter Image: Construction of the second s	Time Filter Image: Control of C	Image: Find state of the s	Image: Finite Filter Image: Filter Ima	Image: Find state of the s	Image: Filter Image: Filter<

Press the green Save button once the desired Date and Time ranges are input.

The Audit Trail screen will then display only messages generated during that Date and Time range.



Settings

Settings are used to create locations, groups, presets, sample plans and certifications. In addition, instrument settings, connectivity, and users are managed from the Settings screens.

The Settings gear button will only be available to "Admin" and "Power User" level User logins. When an "Operator" level username is logged in, the Settings gear button will be disabled and the operator logged in will not have access to any of the settings screens. If Users are disabled then the Settings button and screens will be active and accessible to any user.

From the Home screen press the Settings gear icon to display the settings screen button list for Location, Preset, Sample Plans, Certifications, Users, Connectivity, Instrument, and Service.

Select a settings button from the right column and that settings screen will highlight blue to indicate it has been selected and is being displayed. In the example below the Location button is highlighted blue and the Location screen is displayed.




Blank Page

Instrument Settings


The Instrument Settings screen allows the user to view and edit options and settings for date/time, audio/visual, global channels, warning/alarm data count format and units, sample and zero count settings, data output and formatting, and saving and loading configurations files.

Access the Settings Screen

From the Home screen press the gear icon to display the settings screen.




Home Screen



Access Instrument Settings

From the Settings screen press the Instrument button to display the Instrument Settings screen.



Group Name Location Name 0 Location Add New Group Create New Locations (+) (+)123 Pres \uparrow Group 1 Location 0001 Group 2 Location 0002 Sample Plan Group 3 Location 0003 Certification Location 0004 Location 0005 . Users Location 0006 Location 0007 Connectivit 🗾 Edit Delete 🗾 Edit Delete 3/250 250 / 2500 . Sen

Settings: Groups and Locations

Configuration

Locatio

•

Time and Date

Press the Time button to edit the date, _ date format, time zone, time, and time format, and to establish NTP (Network Time Protocol).

Network Time Protocol

The ApexZ supports a network protocol for clock synchronization called NTP, Network Time Protocol. This service is able to synchronize the time subsystem in the unit based off a NTP server input by the user.

 \uparrow Date 123 Prese 🚞 Date 05/04/2021 Sample Plan Date Format MM/dd/yyy Certificatio Time User Network Time Protocol (NTP) Server Connectivit \$ Instrumer 8 Servi

Instrument Settings Screen

Output

Options

The ApexZ will then query that server for the correct date and time and maintain itself accordingly.

Time

Enable the use of NTP by checking the checkbox in the Time tab.

Press in the empty box to bring up the Enter NTP Server keyboard. Enter a viable NTP server name and press Apply.

Audio / Visual

Set the NTP Poll interval to a value between 64 and 9999. This will be how often in seconds that the unit will query the input NTP server in order to maintain proper timekeeping.

As long as the NTP option is checked, the Date and Time fields in the Time screen will be non-interactable. The user will still be able to change the Date Format, Time Zone, and Time Format.

Set Current Date

Note: If the checkbox on the Network Time Protocol line is checked, the user will be unable to select the Date button.

Pick Date (month/day)

Press the Date button to display the Pick Date screen.

Select today's day, month, and year.

Select the green Save button to save the date settings.




Or press the red X button to cancel.





Set Date Format





Press the green save button to save the date format setting.

248083466-1 R11

Set Time Zone

Scroll down in the Time screen and press the Time Zone button to display the Select Time Zone screen.

Time	Audio / Visual	Options	Output	Configuration	1	•	Location
LONGWOIK	ine Protocor (NTP)			· · · ·	\uparrow		
4.≣	NTP Server		time.goo	ogle.com		_	Preset Sample Plans
4 6	NTP Poll Interval (se	:0)	Í	64		-	Certifications
Time Zone	•	Am	erica/Los_	Angeles		×	Users
Time			3:0)2:47 pm		윪	Connectivity
CTime Form	at			12 HR	\mathbf{V}	۵	Instrument
☆ ≡	: 🔒 🌣					2	Service

Instrument Settings: Time Screen

Press and select your UTC-Time Zone.

Press the green save button to save the
UTC-Time Zone.



Select Time Zone screen

>	<		Select	Time Zone			
		Time	Zone	Location		1	Click and drag the scroll bar
->	υтс	-08:00 ((PST)	America/Los_Angeles	\uparrow	1	
	υтс	-08:00	(PST)	America/Vancouver		1	to scroll up and
	υтс	-08:00 ((PST)	America/Tijuana			down the Time
	UTC	-07:00 ((MST)	America/Denver			Zone List.
	υтс	-07:00 ((MST)	America/Boise			
	υтс	-07:00	(MST)	America/Phoenix			Or click on the
	υтс	-07:00 ((MST)	America/Mazatlan			arrow to scroll.
	υтс	-07:00 ((MST)	erica/Dawson_Creek			
	UTC	-07:00 ((MST)	America/Hermosillo	\mathbf{V}		

Set Current Time

Instrument Settings: Date and Time Screen


Note: If the checkbox on the Network Time Protocol line is checked, the user will be unable to select the Time button.

Press the Time input button to display the Enter Time keypad.

Enter the current time from the numeric keypad as HH:MM:SS Select am or pm from the toggle button.

Press the green Apply button to save the current time entered.

Apply




Enter Time keypad


×	1:55:2	20 pm	
Enter Time		(pm	00:00:00
1	2	3	delete
4	5	6	
7	8	9	
			Apply

Set Time Format

Press the Time Format button to toggle between 12 or 24 hour time format.



Instrument Settings: Time Format



Audio/Visual

Volume

Touch and slide the volume button from 0 - 100%.

Brightness

Language

-

Select the desired language.

Touch and slide the brightness button 0 - 100%.

Press the green save button.

Note: Enabled Alarms will always sound at 100%.

Press the Language button to display the Select Language screen.

Warning: Changing the instrument language will update the instrument firmware and the printer. It will take 1-2 minutes for the update to complete. During this time make sure the ApexZ instrument has at least 50% battery life remaining or connected to the AC adapter and

do not turn off the instrument while updating. Failure to follow these instructions could

Instrument Audio/Visual Screen

Options

Audio

Time

Volum

Audio / Visual

Output

Configuration

 \uparrow

npie Pla

Visual Brightness Uisual Certifications Connectivity Con

Select Language Screen



Decimal Mark

Press the Decimal Mark button to display the Select Decimal Mark screen. Select period or comma. Then press the green save button.

cause the printer to become damaged and unusable.

Digit Grouping

Press the Digit Grouping button to display the Select Digit Separator screen. Select none, space, comma or period. Then press the green save button.

Select Decimal Mark Screen



Select Digital Grouping Screen




Options

Lock Global Channels


ApexZ can lock specific particle size channels globally for all manual samples, preset samples, and sample plans. For example the user can lock channels globally for the 0.5µm and 5.0µm channels and then for all new samples (manual, sample plan, or presets) only 0.5µm and 5.0µm channels will be sampled and those particle size channels counted. Certification Samples are not affected by Locked Global Channels being enabled. Certifications will always use the channel sizes required from the certification standard.

How to Lock Global Channels

Click on the Global Channels buttons to display the edit Global Channels pop up screen. _



Global channels highlighted in green are ENABLED. Global channels highlighted in grey with strike through text are DISABLED. Press any green channel button and it will highlight grey with strike through text to indicate this channel has been disabled. Press any grey button and it will highlight green indicating it has been enabled.




Global Channels Screen

Note: If Lock Global Channels is <u>disabled</u> then all particle channels will be counted and recorded.



Global Channels



Select either Raw or Norm (Normalized) particle counts for the Warning and Alarm limits.

Raw displays total particles counted over time.

Norm displays normalized counts per unit of measure. -

If normalized is selected then select normalized units: $ft^3/m^3/L$.

Tool Tips

Click on the icon/label "Lock Global Channels" and a brief description of that feature will appear.

Tool tips are included for icon/labels included in: Instrument: Options, Output and Configuration Connectivity: AD/LDAP and CIFS input screens



Tool Tips: Lock Global Channels



User Options

Enable Users (enable/disable)

Press the Enable Users button and go directly to the Login Screen and require a User Login. Note to login to ApexZ at least one "local" Administrator level user must have been created prior.

AD Groups (enable/disable)

Enable AD Groups to assign active directory groups to an ApexZ permission levels.

Local Password Cache (enable/disable)

Enable local password cache to save passwords locally on ApexZ when Active Directory Server is unavailable.

Enable Guest Account (enable/disable)

Press to enable guest account and then a guest user will be permitted to login without a password and will be assigned the default "Operator" level read only access. Only an Admin level user can enable the guest account feature.

Auto Lock (enable/disable)

Press to enable Auto Lock and automatically lock the screen after the timeout period.

Timeout

Enter the timeout period to lock the screen when auto lock is enabled.

Include Sample Time (enable/disable)

If Include Sample Time is enabled then the timeout timer will countdown during idle screen time or while sampling. If "Include Sample Time" is disabled then the timeout timer will not count down during a sampling session.

Login Retry (enter login retry limit)

Set the number of times ApexZ will allow a user to retry failed login attempts. After the login retry limit has been exceeded ApexZ will lock out further login attempts and start the Account Lock Duration countdown time. During the Account Lock Duration time period ApexZ will not allow user login. After the Account Lock Duration time period has been exceeded then ApexZ will allow user login again.

Account Lock Duration (enter account lock time duration period)

After login retry is exceeded, ApexZ user login will be locked for the Account Lock Duration time period set here.



Instrument Options Enable Users

	ь 2 0р		ock Duration	00:30:00	\downarrow	۵	Instrument
↑ ≡	-	٠				s.	Service

Home Screen - Auto Locked

Unlock screen from Auto Lock

Press the User Button to unlock user screen.

Then select

- Switch User (new login required)
- Unlock (password required)
- Cancel (remain locked)



	ны	ne s	CIE	en.	·Au		UCKEU	
■ 05/22/2019 10:12:44 am	Admin 1	ψ		ۍ ۲	0ŧ	ů	() •	Group 1 Location 0001
μm 0.5				Σ 0			≡,	
5.0				0			E .	Sample Plans
							Mode: Cycles: Delay: Sample: Hold: Recs:	Manual Sample 0 / 1 00:00:00 00:01:00 00:00:00 9 / 10000
								IDLE
↑ ≡								1
					l.	اء ما	al	

Note: All buttons disabled during auto lock.


Auto Pump Start

Sample Options

With Pump Startup Disabled: The ApexZ pump will start up after the sample Delay time has expired. With Pump Startup Enabled: The ApexZ pump will start immediately after pressing the sample Start button.

Disable Sample Abort

When Enabled, "Disable Sample Abort" will not allow a sample that has been started to be stopped. When Disabled, any sample may be stopped at any time.



Zero Count Options


Mode: Press the Zero Count Mode button to select one of the following 3 modes.

Disabled: The zero count button will be inactive and will not display.

Manual: Zero count button will display and run for the sample time when pressed.

Auto: In Auto mode, zero counting will be successful after sustaining a full minute of no counts. During the test, the timer progress bar will reset to 0 if a count is detected.

Zero Count Mode





Then press the green save button.

248083466-1 R11

Instrument Options: Sample

Test Duration

Press the Test Duration button to enter your desired zero count sample time with the pop up numeric keypad.

Delay Time

Press the delay button to enter the desired delay time.

Delay Time will start after pressing the zero count start button on the Home screen. The Zero Count Sample time will start immediately after the delay time is completed.

		Time	Audio / Visual	Options	Output	Configuration	•	Location				
	•	Channel		123	Preset							
ľ	••	User	Sauto Pump Si	Auto Pump Start								
	Ĩ		💐 Disable Samp	CDisable Sample Abort								
	8	Sample		Zero Count								
	E I	Data Tags	EMode			Manual	<mark>.</mark>	Users				
	=	Audit Trail	• Test Duration			00:00:06	윪	Connectivity				
		omments	🕙 Delay			00:00:00						
							\$	Instrument				
R.	ŵ	=					2	Service				

Instrument Options: Zero Count

Data Tag Name Options

Data Tags are set up in the Instrument Settings \rightarrow Options \rightarrow Data Tags screen. ApexZ supports up to 5 unique Data Tags. **NOTE:** Adding a new Data Tag will require clearing all data records.

Data Tag Name

Any Data Tags that are set up will be listed here by the name given during setup.

Prompt

Each listed Data Tag will also display whether it is set to trigger Pre- or Postsample

Annotations

Each listed Data Tag will show how many Annotations are set up for that Data Tag.

Add New Data Tag

Press the Add New Data Tag button and , the New Data Tag screen will display.

Prompt: Enable the Prompt feature to force a prompt to appear and request to enter a data tag value.

When: Select Either Pre Sample or Post Sample.

Select "Pre Sample" and the data tag _____ prompt will appear when the start button is pressed before a sample starts.

Select "Post Sample" and the data tag prompt will appear after the sample has completed and stopped.

0 Locatio Audio / Visual Configuration Options Output Time 123 Channe Pres Data Tag Name Prompt Annotations Add New Data Tag Sample Plans Use Pre Sample test tag Batch1 Pre Sample Certification: Sample post tag Post Sample Data Tags User Audit Trail Connectivit 🗾 Edit Delete * Instrumen Data Tags: 3 / 5 2 Service



Allow Bypass: Enable this feature to allow the user to bypass entering the data tag prompt when it appears. User will be allowed to enter the data tag value later in the data record screen. *Note: Bypassing a data tag prompt will be noted in the Audit Trail.*


Apply to Certifications: Enable Apply to Certifications to prompt for data tag entry when running certifications.

Lock Entry: Enable Lock Entry to lock the data tag entry on the data record the first time it is set. With Lock Entry enabled no data tags may be modified once they have been set.

Note: If a data tag has been bypassed, user can still set the data tag value in the data record screen. Once the data tag has been set, then the value is locked and cannot be modified.

Barcode Scanner

ApexZ supports data tag input through a standard USB barcode reader connected through the ApexZ USB port. Once the USB scanner is connected, a barcode icon will appear on the home screen header.



Data Tag Annotations

Up to five unique Data Tag Annotations can be New Data Tag × set and used with each Data Tags. Option Saved Data Tag Annotations Press the Data Tag Annotations button to Add New Annotation Data Tag notation create new Data Tag Annotations. User 4 Nam Press the Add New Annotation button to bring up the Enter Data Tag Annotation keyboard and name the Annotation. 🗾 Edit Delete Enter Data Tag Annotation w Shift Spac Apply

Lighthouse Worldwide Solutions

Data Tag Name	×					Nev	/ Dat	la T i	ag	_			
Press the Name tab and	Enter Tag Nam	e									Batch	า# 🧧	Options
enter the new Data Tag name using the keypad /	E.	1			#	%	,	^		+	Enter Valid I	lame 😑	Data Tag Annotations
			<u>\</u>	1	~	<	>	2	3	•		-	Name
	123			Γ	?	Γ		•	•		Delete		
	ABC				s	pace				Ap	riy		
Then press the green Apply button.													

Enter New Data Tag Name Screen

The screen will change to the New Data Tag confirmation screen.


Press the green Save button to save the new data tag name.


NOTE: Adding a new Data Tag will require the user to delete all records currently on the unit.

×						Ba	atch#						
Enter T	īag Nam	e											Options
]	1			#	%	^	*	+		, r	-	Data Tag Annotations
		an a	N	T	۲	< >	> 2	а	٥			•	Name
	123		•		?	1	11	•		Delete			
	ABC				Sp	ace			Ap	ply			


Save New Data Tag Name Screen

ApexZ Operators Manual






Delete Data Tag



Warning: Data records will be cleared when any data tag is added, data tag name is edited, or deleted while there are data records in the data buffer. Press the red X button to cancel and not delete any data or press the green check button to confirm, and data records will be cleared to continue.

Audit Trail Comments

The Audit Trail Comments screen allows the user to pre-define comments for use in the Audit Trail. Comments set up here will be listed in the "Saved Audit Trail Comments" list when adding a comment to a message in the Audit Trail. The user will be able to select these saved comments in lieu of typing in the comment every time.





Add a New Comment

Press the Add New Comment button to bring up the Enter Audit Trail Comment screen.

Enter the text desired for the comment _____ and press the green Apply button to save the comment.

Each saved comment will display on the page. To Edit or Delete them, first select the comment then press the Edit or Delete button.





Output

The Output screen allows the user to adjust output settings and data formatting for USB, Print, and CIFS.

	Instrument Options: USB Output
USB	Time Audio / Visual Options Output Configuration Classifier
Press the USB button to display the Data Output and	Print File Format
Format options for USB.	Certifications
Data Record File Format	
Set the file format to .csv, .xlsx, .lsdx, or .pdf.	
Data and records output to USB will be saved in the	CFile Format
selected format.	instrument
On Sample (.csv format only)	🛧 🔳 🖹 🗢 🔨 Service
Enable Output "On Sample" and output will initiate after	avany sample

Enable Output "On Sample" and output will initiate after every sample.

On Alarm (.csv format only)

Enable Output "On Alarm" and output will initiate after every alarm is exceeded.

Certification Report File Format

Set the file format to .lsrx or .pdf. The lsrx file format is a non-human readable format that can be imported into LMS XChange and contains both the Certification Report and the data records used to generate the report.

Print – On-Board Printer

The ApexZ can print to either the optional embedded front-panel printer or to a network printer. The first screen in the Print tab covers automatic print options for the optional on-board printer.

On Sample

Enable Output "On Sample" and output will print – after every sample is completed.

Average

Enable Average to print the average counts of all cycles at a given location after all cycles have completed.

On Alarm

Enable Output "On Alarm" and output will print / after every alarm is exceeded.

Instrument Options Output

/	Tim	18	Audio / Visual	Options	Output	Configuration	•	Location
	÷	USB	On Boo	ard Printer	Net	twork Printer	123	Preset
	÷	Print	100				=,	Sample Plans
	12	CIFS	h ‡ Ave لمن	rage			e ,	Certifications
	For	matting					*	Users
	/						器	Connectivity
							*	Instrument
	â		a				2	Service


Print – Network Printer

The ApexZ is able to send print data to a properly set up network printer.

Select the Network Printer tab from the Print screen.

Record Format

The Network Printer output can be set to either Printer Emulation or Data Table



Printer Emulation will display the data in three columns on each page, retaining the size and format of what would come out of the onboard printer.

Data Table will display the output information in a table format.



CIFS

Common Internet File System

Select the CIFS button.

Data Record File Format

Select the .csv, .xlsx, .lsdx (Lighthouse Secure Data Transfer), or .pdf file format.

Include Audit Trail

Select Include Audit Trail and all audit trail records will be included in the output.

	Ti	ne	Audio / Visual	Options	Output	Confiç	juration	•	Location
	÷	USB		123	Preset				
	÷	Print	File Format .Csv .xlsx .lsdx .pdf						Sample Plans
•	5	CIFS		eş.	Certifications				
		rmatting		🖹 Certif	ication Rep	oort		*	Users
			🔒 File Format			.lsrx	.pdf	윪	Connectivity
								*	Instrument
	A							2	Service

Instrument Output CIFS

Clear Data Buffer

Select Clear Data Buffer (only valid when exporting ALL data) and after the data file is successfully copied to the CIFS share the ApexZ instrument database will be cleared. If the data exporting fails then the ApexZ database will not be cleared.

Certification File Format

Select .lsrx or .pdf and Certification outputs to CIFS will be saved in the highlighted format. The lsrx file format is a non-human-readable format that can be imported into LMS XChange and contains both the Certification Report and the data records used to generate the report

Output Formatting - Print

Print Output Format -

Press the Print Tab to edit formatting parameters for printed records and certifications.

Cumu/Diff

Select to print data in either Cumulative, Differential, or both.

Raw/Norm

Select to print Raw or Normalized data. When selecting Normalized also select the Normalized Units (cubic feet, cubic meters, or Liters).

Reverse

Select Reverse and printing will print from the bottom up.

Font Size

Select Large or Small font size for print output.

--Use the scroll bar to view more formatting options--



File Header

Press the File Header Edit button to enable/disable file headers.

By default all File Header System names will be listed on the left column and highlighted blue indicating they are all selected. A matching Name will appear on the right side of the screen.

Select a highlighted System Name and the highlight will disappear indicating that system name is inactive and the matching Name will be removed and will not appear on the list of Names. Select the inactive System Name and it will highlight blue and a matching Name will appear on the Name list on the right side of the screen.

Adjust Name order on the list

To adjust the order of Names press and hold the dotted box to the left of each name on the list. Drag the name up or down to change its position.

Time	Audio / Visual Options		Output	Configuration	on	•	Location
usb 🕴	Print		.csv / .xlsx	.pdf		123	Preset
📑 Print	File Header		Þ	Edit	\uparrow	I,	Sample Plans
	B+Column Form	mat	Edit				
CIFS	tatistics		E	Edit		E.	Certifications
Formatting	Company					×	Users
		 				윪	Connectivity
	• Address		 		\checkmark		Instrument
						~	Service
			•				

Column : 8	selected : 8
System Name	Custom Name
Company	iiii Company
Address	iiii Address
Model	iiii Model
Flow	III Flow
Serial #	iiii Serial #
Cal Date	IIII Cal Date
Instrument Name	iiii Instrument Name
Separator	iiii Separator

Column Format

Press the Column Format Edit button to enable/disable the columns included in printed data.


By default all System Names will be listed on the left column and highlighted blue indicating they are all selected. A matching Name will appear on the right side of the screen.


Select a highlighted System Name and the highlight will disappear indicating that system name is inactive and the matching Name will be removed and will not appear on the list of Names. Select the inactive System Name and it will highlight blue and a matching Name will appear on the Name list on the right side of the screen.

Adjust Name order on the list

To adjust the order of Names press and hold the dotted box to the left of each Name on the right-hand list. Drag the Name either up or down to change its position.

Any Data Tags that have been created will be listed at the end of the System Name list. These can be included/excluded the same as any of the System Names, and can also be re-ordered in the right-hand column the same as the other Names.





Statistics

Press the Statistics Edit button to enable/disable the columns included in summary Statistics.

By default all System Names will be listed on the left column and highlighted blue indicating they are all selected. A matching Name will appear on the right side of the screen.

Select a highlighted System Name and the highlight will disappear indicating that system name is inactive and the matching Name will be removed and will not appear on the list of Names. Select the inactive System Name and it will highlight blue and a matching Name will appear on the Name list on the right side of the screen.

Adjust Name order on the list

To adjust the order of Names press and hold the dotted box to the left of each Name on the right-hand list. Drag the Name either up or down to change its position.

Company Name


The user is able to input custom text for the Company Name column that can be included in the File Header.


Press in the box under Company Name to bring up the Enter Company keyboard. — Enter the desired Company Name and press the green Apply button.

Address

The user is able to input custom text for the Address column that can be included in the File Header.

Locatio Time Audio / Visual Options Output Configuration Prese .csv / .xlsx Print .pdf File Header Edit \uparrow 🖡 Sample Plans Edit E+Column Format Certifications CIFS Edit 🔒 Statistics Company 2 Users Formatting Connectivity • Address * Instrument J. 8





Press in the box under Address to bring up the Enter Address keyboard. Enter the desired Address and press the green Apply button. Up to 128 characters can be used in the Address field.

Neither the Company Name nor Address support multiline entry.

Output Formatting - .csv/.xlsx

.csv / .xlsx Output Format

Press the .csv/xlsx tab to edit formatting parameters for the .csv and .xlsx file types.

Raw/Norm

Select to print Raw or Normalized data. When selecting Normalized also select the Normalized Units (cubic feet, cubic meters, or Liters).

File Header

Press the File Header Edit button to enable/disable or edit custom name data for the file headers in data output to .csv and .xlsx files.

By default File Header System Names will be listed on the left column with names highlighted blue indicating those are selected. A matching Custom Name will appear on the right side of the screen.

Some System Names may be grey and not highlighted blue. This indicates they are available but not selected for output.

Select a System name to enable (highlighted blue) or disable (grey no highlight). Enabled System Names will appear on the Custom Name list on the right side of the screen. Disabled System Names will not appear on the Custom Name list.

Adjust Custom Name order on the list

To adjust the order of Column Format Custom Names press and hold the dotted box to the left of each custom name on the list. Drag the custom name either up or down to change its position.

Edit Custom Name

Press the Edit button to display the Enter

Locati 0 Audio / Visual Options Configuration Time Output 123 USB Pres Print .csv / .xlsx .pdf 🔒 Raw / Norm Print Sample Plan Raw Norm Certification CIFS Sormalize Units mª 8 User Formatting Edit File Header Column Format Edit Connectivity .xlsx File Password * Instrumen 8 Servic

×				File	Head	er					•	
	Colu	mn : 7						select	ed:7			
System Name			Custom Name									
Created:			IIII Created:									
Created By:						Crea	ted By:					
Model Name:						Mode	el Name:					
Serial Number:					iiii Serial Number:							
Flow Rate:						Flow	Rate:					
Flow Units:						Flow	Units:					
Instrument Name:						IIII Instrument Name:						
							<u>,</u>		Edit			
×				Flo	w Rat	e:						
Enter Custom Nam	ne								F	low F	Rate	
q	w	e	r	t	у		u		o	р		
	а	s	d	f	g	h	j	k	1			
Shift	z	×	с	v	Ь		n	m		Delete		
123					Space				Apply			

Custom Name keypad. Enter a new custom name and press the green Apply button to save the custom name to the Custom Name list.

Column Format

Press the Column Format Edit button to enable/disable or edit custom name data columns output to .CSV and .xlsx files.

By default Column Format System Names will be listed on the left column with names highlighted blue indicating those are selected. A matching Custom Name will appear on the right side of the screen.

Some System Names may be grey and not highlighted blue. This indicates they are available but not selected for output.

Select a System name to enable (highlighted blue) or disable (grey no highlight). Enabled System Names will appear on the Custom Name list on the right side of the screen. Disabled System Names will not appear on the Custom Name list.


Adjust Custom Name order on the list

To adjust the order of Column Format Custom Names press and hold the dotted box to the left of each custom name on the list. Drag the custom name either up or down to change its position.

Edit Custom Name

Press the Edit button to display the Enter Custom Name keypad. Enter a new column format custom name and press the green Apply button to save the custom name to the column format settings.





× Sample (sec) *											
Enter Cu	stom Name								Sar	nple (sec)
	q	w	e	r	t	у	u	i		р	
		a	s	d	f	g h	j	k	1		
	Shift	z	×	с	v	ь	n	m		Delete	
	123				Space			Apply		ply	

.xlsx File Password

The ApexZ enables greater security for .xlsx (Excel) files via the addition of an optional file password. This function will add a password to all Excel files created, prompting any user to enter that password before attempting to edit the file.



Press the gray box on the File Password line to bring up the File Password screen.

Use the keyboard to enter the desired password, and verify the password by entering it a second time. Press the green Apply button to verify the entries.





Output Formatting - .pdf

Press the .pdf tab to configure the parameters for data output as a PDF.

There are two layout options: Printer Emulation and Data Table.

Printer Emulation will lay out the included data in columns emulating what would come out of the optional on-board printer. An example is shown in the screen.





Regardless of the layout chosen, the formatting options remain the same and will affect both layouts.

Raw/Norm


screen.

options--

Select to output Raw or Normalized data. When selecting Normalized also select the Normalized Units (cubic feet, cubic meters, or Liters).

Cell Align

Choose whether to align the text in the cells to the Left, Center, or Right.



File Header

Press the File Header Edit button to enable/disable or edit custom name data for the file headers in data output to .CSV and .xlsx files.

By default File Header System Names will be listed on the left column with names highlighted blue indicating those are selected. A matching Custom Name will appear on the right side of the screen.

Some System Names may be grey and not highlighted blue. This indicates they are available but not selected for output.


Select a System name to enable (highlighted blue) or disable (grey no highlight). Enabled System Names will appear on the Custom Name list on the right side of the screen. Disabled System Names will not appear on the Custom Name list.

Adjust Custom Name order on the list


To adjust the order of File Header Custom Names press and hold the dotted box to the left of each custom name on the list. Drag the custom name either up or down to change its position.

Edit Custom Name

Press the Edit button to display the Enter Custom Name keypad. Enter a new custom name and press the green Apply button to save the custom name to the Custom Name list.







Column Format

Press the Edit Column Format button to enable/disable or edit custom name data columns output to PDF files.

By default Column Format System Names will be listed on the left column with names highlighted blue indicating those are selected. A matching Custom Name will appear on the right side of the screen.

Some System Names may be grey and not highlighted blue. This indicates they are available but not selected for output.


Select a System name to enable (highlighted blue) or disable (grey no highlight). Enabled System Names will appear on the Custom Name list on the right side of the screen. Disabled System Names will not appear on the Custom Name list.


Adjust Custom Name order on the list


To adjust the order of Column Format Custom Names press and hold the dotted box to the left of each custom name on the list. Drag the custom name either up or down to change its position.

Edit Custom Name

Press the Edit button to display the Enter Custom Name keypad. Enter a new column format custom name and press the green Apply button to save the custom name to the column format settings.









PDF File Password

The ApexZ enables greater security for PDF files via the addition of an optional file password. This function will add a password to all PDF files created, prompting any user to enter that password when trying to open the file in a PDF editor.

WARNING: Once the PDF file password is set and applied, it cannot be removed without an instrument settings reset or factory reset. If the password is added be sure to record it well.




Use the keyboard to enter the desired password, and verify the password by entering it a second time. Press the green Apply button to verify the entries.



Press the gray box on the File Password line to bring up the File Password screen.

Configuration

All of the ApexZ settings are saved in the ApexZ configuration file and can be exported onto a USB key. The exported configuration file can then be loaded onto another ApexZ instrument or edited using our LMS XChange software.



Save ApexZ configuration file to USB key

Place a USB key in the ApexZ USB drive. Press the Save button to display the Save Confirmation screen.


Instrument Options: Configuration

Time	Audio / Vis	ual Optic	ins	Output	Configu	ration	Location
		Cont	figuration			15:	B Preset
Current						Save	Sample Plans
Load Co						Load	Certifications
							Users
						몲	5 Connectivity
						8	Instrument
^ :		٠				4	Service


Press the green check mark button to confirm exporting the configuration to the USB key.

Or press the red X button to cancel.

Save Configuration File



Load ApexZ configuration file from a USB key



Blank Page

Connectivity

Special Note:

LMS XChange Software 2.x

The LMS XChange 2.x software was designed for use with ApexZ particle counters and is downloadable from the <u>www.golighthouse.com</u> website. LMS XChange 2.x allows users to import data collected on ApexZ3 and ApexZ50 portable particle counters and display the data in a table that the user can view and/or save to an Excel, .CSV or PDF file. Additionally, Administrator users can manage configurations for ApexZ instruments using LMS XChange. LMS XChange has the following main features:

- Import data manually from Lighthouse ApexZ3 and ApexZ50 particle counters.
- Import .lsdx (Lighthouse secure data transfer) files created by ApexZ3 and ApexZ50 instruments.
- Display the imported data in a table, allowing the user to also export the data to an Excel, .CSV or PDF file.
- Three levels of user accounts Operators, Power Users, and Administrative users.
- LMS XChange Administrator users can create, edit, save and load ApexZ3 and ApexZ50 configurations.
- Record logins, updates, delete, import, configuration changes, etc. in the LMS XChange Audit Log.


From the Settings screen press the Connectivity button to view Ethernet, WiFi, AD/LDAP, CIFS, Network Printer, and 3rd Party API settings.




Connectivity


Network - Ethernet

Connectivity: Ethernet Enabled



Connectivity: Ethernet Disabled





Connectivity: Ethernet

IP Address



Enter the valid Address and press the Apply button.

When IP, Subnet and Gateway are entered then press the Enable button to enable the Ethernet connection

Subnet Address

× 255.255.255.255										
Enter Subnet										
1	2	3	delete							
4	5	6								
7	8	9								
	0		Apply							

Gateway Address

× 10.9.8.7								
Enter Gateway								
1	2	3	delete					
4	5	6						
7	8	9						
	0	•	Apply					

DHCP

(Dynamic Host Configuration Protocol)

When DHCP is enabled the IP Address, Subnet, and Gateway will be dynamically assigned from the HOST network connection.

User editing will be locked.

To enable DHCP, disable the Ethernet connection first.

After the DHCP has been enabled, press the Enable Ethernet button to complete the DHCP connection.

Network	AD / LDAP	CIFS	Network Printer	3rd Party API	•	Location
e Ethernet	T Enable				123	Preset
	IP Address			5.5.5.5	II. s	ample Plans
🔶 WiFi	& Subnet			6.6.6.6		
	Gateway			7.7.7.7	.	Certifications
	MAC Address		88:88:88:88	8:88:88	*	Users
	рнср				윪	Connectivity
					۵	Instrument
^					2	Service

Connectivity: DHCP enabled

Connectivity: Wi-Fi Screen

Network - Wi-Fi


The IP Address, Subnet, and Gateway for Wi-Fi can be either manually entered (BEFORE enabling the Wifi) or assigned via the DHCP.

Edit Wifi Address

With the Wifi and DHCP checkboxes disabled, press the IP address, subnet or gateway input buttons to display the input key pad and edit these values.

Network	AD / LDAP CIFS	Network Printer	3rd Party API	•	Location
e Ethernel	? Enable	>	<u> </u>	123	Preset
- WiFi	P Address	10.8.125	.100		Sample Plans
₩IFI	Subnet	255.255.1	28.0		Sample Flans
	🕎 Galeway	10.8	3.0.1	.	Certifications
	MAC Address	44:44:44:44	44:4	*	Users
	量онср			윪	Connectivity
	Priority SSID	Sta	tus 🗸	*	Instrument
	1 (+) Ad	d New Access Point	*		
🔶 🖿				~	Service

IP Address



Subnet Address



Gateway Address

Enter Gateway								
1	2	3	delete					
4	5	6						
7	8	9						
	0		Apply					

Enter the valid Address and press the Apply button.

When IP, Subnet, and Gateway are entered, press the Enable button to enable the Wifi connection

248083466-1 R11

DHCP

(Dynamic Host Configuration Protocol)

When DHCP is enabled the IP Address, Subnet, and Gateway will be dynamically assigned.

User editing will be locked.

To enable DHCP, disable the Wifi connection first.

After the DHCP has been enabled, press the Enable Wifi button to complete the DHCP setup.


Network	AD / LDAP	CIFS	Network Printer	3rd Party API	•	Locat
u Ethernet	Enable			\uparrow	123	Pre
	P Address	2	255.255.255	.255	=	Sample Pl
<u>~</u>	Subnet		6.	6.6.6		
	Gateway		7.	7.7.7	-	Certificati
	MAC Address	4	14:44:44:44	:44:4	*	Us
					윪	Connecti
Р	riority Priority SSI	D	Sta	itus	\$	Instrum
	<u>^</u>	🕀 Add Ne	w Access Point	\downarrow		mətrum

Connectivity: DHCP enabled

Add New Wifi Access Point

With the Wifi disabled, the user is unable to add new Wifi networks to the connection list.


Enable Wifi in order for the "Add New Access Point" button to be enabled.



After a Wi-Fi has been selected, the ApexZ will mark with a red dot the remaining required tabs to complete the configuration.

Open Wi-Fi configurations: Security will be preset to "None". Local Network Wi-Fi: Security will be preset to WPA-PSK (also supports WPA2). Password input is required.

Enterprise Wi-Fi configurations: Security will be preset to WPA-EAP (also supports WPA2). Username input is required.



Select Wi-Fi Security

Security will be selectable only for Manual Wi-Fi configurations.

Security will be preselected according to the Wi-Fi type detected.

WPA-PSK (also supports WPA2) Wi-Fi Protected Access – Pre-Shared Key Password required only.

WPA-EAP (also supports WPA2) Wi-Fi Protected Access – Enterprise Certificate accepted, not required. Username required. Password required.

Connectivity: Wi-Fi Security



Enter Username

Enter Username with the key pad then press Apply.

Then press Apply.



Press the Password Button.





Enter Password

Enter password on the keyboard screen and press Apply.




Press the Save button to save settings and connect to the selected Wi-Fi.





Connectivity: Wi-Fi Password



When the ApexZ successfully connects to the Wi-Fi network, a "Connection Successful" message will appear.

×

All configured Wi-Fi networks will be listed in the Network → WiFi screen, and when selected a Connect box will be shown (if not already connected). This allows users to connect to previously configured WiFi networks without setting them up again.

A Wi-Fi network to which the unit is connected will display in the Status column.




Active Directory AD/LDAP

Validate Users through your network Active Directory

ApexZ allows user login authentication to happen through your Active Directory network connection. User names and passwords logging into the ApexZ instrument will be verified against Active Directory usernames and passwords.

User Permission Levels can be assigned to AD Groups on your network. Users logging in from AD Groups will be given that AD Groups' assigned permission Level on the ApexZ instrument.

Contact your IT department to obtain the URL and Domain for your Active Directory.



Active Directory URL (optional)

Press the URL button and use the keypad to enter the Active Directory URL.

Note: If URL is entered, the unit only attempts to authenticate to the Active Directory server at that IP.

Active Directory Domain (required) Press the Domain button and use the keypad to enter the Domain name.

Note: If Domain is entered and IP is left blank, then the unit will pull all valid IP addresses of the domain and attempt to validate credentials until one succeeds.

Secure

Press to enable the Secure button to use LDAP over SSL protocol.

Test

Press the Test button to test login to the Active Directory.


Note: Active Directory URL and Domain must be set up and a Test Login performed before logging in with AD user accounts.



Connectivity: Active Directory


To activate AD Groups go to:

- 1. Settings
- 2. Instrument
- 3. Options
- 4. User
- 5. AD Groups



CIFS -Common Internet File System

The CIFS "Common Internet File System" is used to connect and export data to shared network drives.



Connectivity: CIFS

Press the Path button and enter the CIFS Path using the keypad.

Then press the green apply button.



Connectivity: CIFS URL



Connectivity: CIFS path



Folders 1-2-3-4

Connectivity: CIFS Folders

Folder 1, Folder 2, Folder 3, and Folder 4 are Network AD / LDAF CIES Network Printer 3rd Party API designed to allow the user to organize data into sub- \uparrow Pre folders according to the instruments they were Sample Pla downloaded from. Folder 1 Folder 2 Certifica Folder 3 2 Hs Folder 4 옮 Connectivi Press each Folder button and select the desired sub Auto Authenticate Test Authentication Y folder name from the pop up list: . Serial # 6 Select Folder 1 Date/Time Model Serial # + Folder 1 Instrument Name Serial # Date/Time None Model Instrument Name None

After selection, press the Save button and the entry will be listed for that Folder. The buttons for subsequent Folders will only be available once the previous Folder has a selection saved.

The Instrument Name folder will be named after what is entered on the Settings->Service->About screen, in the Instrument Name field. If this field is empty, selecting Instrument Name for a CIFS Folder will show an orange folder icon in the field after saving, indicating that the Instrument Name has not yet been set on the Apex unit.



Auto Authenticate

Press the enable Auto Authenticate button and ApexZ will use the currently logged-in user's credentials to authenticate with the network drive when exporting data. Disabling Auto Authenticate will require every user to enter their username and password prior to exporting data to the network drive.




Test

Test Authentication

Press the Test button to test Authentication.

Network Printer



Press the Test button to test the connection to the Network Printer.

3rd Party API - Application Programing Interface

Note: Contact Lighthouse for documentation about API commands to integrate software with the ApexZ

Port Number

Press the Port button to enter a port number with a minimum of 4000 and a maximum of 8000.

×	40	00	
Minimum: 4000 Maximum: 8000			
Enter Port			
1	2	3	delete
4	5	6	
7	8	9	
	o		Apply
	SI ()		

Connectivity: API Port Number



Connectivity: Enable 3rd Party API

Enable API

After selecting the port number then press the Enable API button.

After the API has been enabled the port number may not be modified. To change the port number you must disable the API enable button first to allow editing of the port number.



Blank Page

Locations and Groups


The ApexZ instrument has the capacity to store up to 2500 location names. The first 250 location names have been named as "Location 0001", "Location 0002" and up to "Location 0250". Additional Location names can be added and the location name will default to the next sequential location # for example "Location 0251". Location names can be easily edited/renamed or deleted.

Groups may also be created and locations may be assigned to Group names. Up to 250 Groups can be created and locations can be assigned to these groups.



Home Screen


Configuration: Locations and Groups Screen



Create New Locations

Press the 💽 button to Create a New Location.

After pressing the 🕑 button the next location number will be created and named "Location 0251".



Configuration: Locations and Groups Screen

Add additional locations before saving

Press the 🕒 button to add one more location to the list of new locations. Add Location 7 Location Name Each new location will be named sequentially "Location 0251, Location Add Location \oplus 0252, Location 0253... etc. Ξ Location 0251 . tion 0252 ation 0253 Press the - Delete button to delete a highlighted location from the list. Press the Edit button to change the name of the highlighted location. Z Edit Press the green Save button to add the list of new locations Delete to the ApexZ database.


Add Locations Screen

•

Edit a Location

Select a location name and it will highlight in blue indicating that it has been selected.

Press the Edit button to edit the selected location name or to assign the location to a Preset.



Locations Screen

Edit Location Name

Delete and type in a new location name and then press the green Apply button.


Edit Location Preset

Press the Preset button to display the Assign Preset to Location screen.


Select the desired preset to apply to this location. The selected Preset will display in the Assigned Preset box. Note that a preset setting may be assigned to multiple locations but a single location may only have one Preset assigned to it.

Press the green Save button to save the location edits.

Edit Location Name Screen




Delete a Location



Add a New Group

Groups may be created and locations may be assigned to groups. Multiple locations may be assigned to a group but a single location may only be assigned to one group. Group names can be customized.



Groups and Locations

Adding a New Group

The New Group screen displays a list of all locations available to be assigned to the new group.

Total locations available. -

Total locations selected. -

Note: ApexZ can store up to 250 group names. At least 1 Location must be assigned to a Group and each location may only be assigned to 1 group.

×		New Group		
Locations : 250		selected : 0	•	Location
Location Name		Assigned Locations	ŏ	
Location 0001	\uparrow			Name
Location 0002				
Location 0003				
Location 0004		No Location selected		
Location 0005				
Location 0006				
Location 0007				
Location 0008	\mathbf{V}			

Assign Locations to a new Group

Select each location to be included in the group. Selected locations will highlight blue indicating they have been assigned to the new group and they will display on the Assigned Locations list to the right.

×		New Group	
Loca	tions : 250	selected : 3	 Location
Location	n Name	Assigned Locations	
Location 0001	个	Location 0001	Name
Location 0002		Location 0002	1
Location 0003		Location 0003	
Location 0004			
Location 0005			
Location 0006			
Location 0007			
Location 0008	\checkmark		

Select the Name button to display the Enter Group Name screen.

Group Name

Use the keypad to enter the Group name.



Press the green Apply button to apply the Group name.

Press the green Save button to save and return to the Location/Group List. The green save button only appears after a name has been entered and the Apply button has been pressed.

Group Name Screen

×						New	Grou	р					
Enter G	Group Na	ame							(Grou	р 1	•	Location
	1	2	3	4	5	6	7	8	9	Enter Vali	d Name	•	Name
		-	/	:	;	() \$	&	@				
	#+=				?	!		•		Delete			
	ABC				Spa	ace			Ар	ply			

Group Selected


When a Group name is selected it will highlight blue and the locations assigned to that group will be displayed in the Location Name list.

The Add Available Locations button will be available to add more locations to the Group.

Group Name Location Name Locatio Ð Add New Group Add Available Locations Ð ... Pres Location 0001 Group 1 Location 0002 Sample Plans Location 0003 Certifica . Connectivit 🗾 Edit Delete Edit Delete \$ Instrum 3/2500 1/250 . Servi

Group Selected – Assigned Locations Displayed

No Group Selected – All Locations Displayed



No Group Selected

When no Group Name is selected then all Locations will be displayed and available to Edit or Delete.

The "Add Available Locations" button will change to a "Create New Locations" button and allow the user to create new locations in addition to Edit or Delete locations.

Locations and Groups Screen

Edit a Group

Select a Group name and it will highlight blue and display all assigned locations.


Press the Edit button below the group name list and the Edit Group screen will display.

Select additional locations to add to the Group or select locations to remove. Selected location names appear highlighted in blue.


Highlighted location names appear included in the Assigned Locations list for this group.

Remove a location from the group by deselecting the highlighted location from the location name list and it will be removed from the Assigned Locations list for this group.

Select the Name tab and edit the Group Name.



Edit Group Locations



Edit Group Name




Edit the Group name with the keypad and press the Apply button.



Press the green save button to save the edit.



Delete a Group




Locations and Groups Screen

Delete Group Name Confirmation

Press the green check button to delete the selected group.

Or press the red X button to cancel and not delete.



Blank Page

Presets

Sample setting parameters can be saved and named as Presets. Presets can then be assigned to locations. Locations linked to Presets will be bound by the sample parameters saved in the Preset. Up to 500 Presets can be added, edited, or deleted.

Add a Preset **Preset Screen** From the Home screen Preset Name Delay Sample Time Hold Cycles Locati • press the Settings button. Add New Preset 123 Preset From the Settings screen press the Preset button. 📕 Sample Plans From the Preset screen No Presets Added Certification: press the Add New Preset button. Users Connectivity Edit Î Delete \$ Instrumen 0 / 500 . Servic **Preset Sample Settings**

Mode	Sample			x Delay	: 23:59:59	8	Sample
Delay	00:00:00	Enter V		0:0	0:00	••	Channel / Alerts
Sample Time	00:01:00	1	2	3	delete	🛡 Fr	vironmental
Hold	00:00:00						/ Alerts
Cycles	1	4	5	6		•	Location
Units	m³	7	8	9		-	Name
Volume	0.028		0		Apply		

248083466-1 R11

Preset Sample Settings

Enter Preset Sample Settings.

Delay Time

Select the Delay button and enter the delay time using the keypad and press the green Apply button.

Sample Time

Select the Sample Time button and enter the sample time using the keypad and press the green Apply button.

Hold Time

Select the Hold button and enter the hold time using the keypad and press the green Apply button.

Cycles

Select the Cycles button and enter the number of Cycles using the keypad and press the green Apply button.

Units

Select the Units button and select the units ft3 (cubic feet), m3 (cubic meters), or L (liters).

Volume

Select the Volume button and enter the sample volume with the keypad and press the green Apply button. Note that changing the sample volume will adjust the sample time accordingly. Vice-versa, changing the Sample time will change the sample volume accordingly.

Press the Channel/Alerts button to enable channels, warnings, and alarms.

New Preset × Mode Sample 8 Sample Max Delay: 23:59:59 Enter Value Delay 00:00:00 Channel Alerts ÷ 00:00:00 Sample Time 00:01:00 Environmenta / Alerts delete Hold 00:00:00 Location Cycles 1 Name Units m³ Volume Apply 0.028

Edit Channel

Disable a Channel

Press one of the green colored channel size buttons to disable that channel and that channel will turn a grey color and strike-out the channel size. This channel will not be included in this preset sample setting.

Note if Global Channels are enabled then channel settings may not be edited.

Preset Channel/Alerts Screen

×	New Preset								
Channel	Warning (ft ^a)	Alarm (ft ^a)	Sample Sample						
0.3	0.0	0.0	Channel Alert						
0.5	0.0	0.0	Environmenta / Alert						
1.0	0.0	0.0	• Locatio						
3.0	0.0	0.0	Nam						
5.0	0.0	0.0							
10.0	0.0	0.0							

Preset Warning and Alarm settings

Enter Warning and Alarm settings

Press a warning button for any channel size and it will display the Enable and Enter Warning Screen.

Press an alarm button for any channel size and it will – display the Enable and Enter Alarm Screen.

×	× New Preset								
\sim	Channel	Warning (ft ^a)	Alarm (ft ^s)	Sample					
	0.3	5.0	10.0	Channel / Alerts					
	0.5	0.0	0.0	Environmental / Alerts					
	1.0	0.0	0.0	Location					
	3.0	0.0	0.0	Name					
	5.0	0.0	0.0	•					
	10.0	0.0	0.0						

Enable & Enter Warning



Enable & Enter Alarm

× o					
Enable	Disable				
		10			
2	3	delete			
5	6				
8	9				
0		Apply			
	2 5 8	2 3 5 6 8 9			

Press the Enable button and it will highlight blue indicating the Warning or Alarm has been enabled.

Enter a particle count value with the keypad and press the green Apply button to save that warning or alarm level to this Preset.

The preset warning button will display orange, and the Alarm preset button will display red when enabled.

Environmental Alerts

Choose Temperature Units


Press to select either Celsius or Fahrenheit. The ApexZ will automatically convert the values within the Alert boxes.

Set Temperature Alert Values

Press the value for the Alarm High to bring up the input keyboard. Choose Enable or Disable for the Alert, enter a value (within the range indicated) and press the Apply button.

× 55						
Minimum: 0	Enable	Disable				
Maximum: 150 Enter Alarm H			55			
1	2	з	delete			
4	5	6				
	8	9				
	0	243	Apply			

× 32				
Er.	nablo	Disable		
Minimum: 0 Maximum: 150				
Enter Warning High			5	
1	2	3	dalate	
4	5	6		
7	8	9	. T	
	0		Apply	



Environmental Alert Settings

×	32					
	Enable	Disable				
Ainimum: 0 Aaximum: 150						
Enter Warning	Low		40			
١	2	3	delete			
4	5	6				
7	8	9				
	0	-	Apply			

×	3	32	
Minimum: 0 Meximum: 154	Enable	Disable	
Enter Alarm	Low		35
1	2	3	delete
	5	6	
7	ų	9	
	0	14	Apply

Temperature Alert Settings

Repeat the action for the Warning High, Warning Low, and Alarm Low values.

Set Relative Humidity Alert Values

Press the value for the Alarm High to bring up the input keyboard. Choose enable or Disable, enter a value (between 0 and 100 percent) and press the Apply button.

Repeat the action for the Warning High, Warning Low, and Alarm Low values.





×		0			
	Enable	Disable			
Minimum: 0 Maximum: 100					
Enter Warning Hig	jh		70		
1	5	а	delete		
4	5	6			
7	8	9			
			Apply		

Environmental Alerts Set

×	× 0					
	Enable	Disable				
Minimum: 0 Maximum: 100		2				
Enter Warning	g Low		40			
1	2	3	delete			
4	5	8				
7	8	9				
	0		Apply			

×			
	Enable	Disable	
Minimum: 0 Maximum: 100			
Enter Alarm L	w		
1	2	з	dalete
4	5	6	
7	8	9	
	0		Apply

Relative Humidity Alert Settings

Assign Locations to the New Preset

A preset can be assigned to multiple locations.

Select location

Select location name(s) from the Location Name list and they will highlight blue and display on the Assigned Locations list.

Multiple locations may be selected and assigned to the same preset at the same time.

To remove a location from the assigned locations list press the location name again and it will unhighlight and be removed from the assigned locations for that Preset.

Press the Name button to display the Preset Name Screen.

Preset Name

Use the keypad to enter the preset name.



Press the green Apply button.

You can press any one of the Sample, Channel/Alerts, Environmental Alerts, Loo

Channel/Alerts, Environmental Alerts, Location, or Name buttons and review any of your Preset settings.



Apply

Space

Preset Name Screen



Press the green Save button to save the Preset.



Assign Locations to New Preset

Nam

Locations linked to Presets

When a preset has been assigned to a location then a linked icon will appear before the location name on the Home screen. The link icon indicates the location is linked to a Preset sample setting.


		Home Screen	
	= 07/19/2019 admin 4:01:21 pm	🜵 🔌 😲 🖬 💎	🣀 Group 1
The preset sample setting will be	🌞 _{µm}	Σ	? 🔗 Location 0001
displayed and locked below.	0.3	0	≡ ■ Sample Plans
	0.5	0	-
	1.0	0	123Preset 1Mode:Sample
	3.0	0	Cycles: 0 / 1 Delay: 00:00:00
	5.0	0	Sample: 00:01:00 Hold: 00:00:00 Recs: 0 / 10000
	10.0	0	
			IDLE
	♠	Feed	



Linked Icon

Presets assigned to locations will display the linked icon on the Groups and Locations configuration screen.

Groups and Locations Screen



Blank Page

button to display the Settings Screen.

Sample Plans

Sample Plans are a collection of Locations with custom sample settings assigned to each location. Recipes (supported in previous versions of ApexZ instruments) were a collection of locations with the same sample settings used at each location. Sample Plans can act as Recipes did, but they can also, in addition, assign different sample settings to each location, either via a Preset or manually.

06/04/2019 10:50:24 am U Admin 1 ŵ 04 X 0 Group 1 From the Home screen press the Settings 0 ${\mathscr S}$ Location 0001 ¢ Σ μm =, 0 0.3 Sample Plans 0.5 0 ≝ 1.0 0 Manual Sample Mode: 3.0 0 Cycles: Delay: 00:00:00 5.0 Sample: Hold: 0 00:01:00 00:00:00 Recs: 0 / 10000 10.0 0 IDLE ♠ \triangleright 101

Home Screen

Settings: Sample Plans Screen

Press the Sample Plans button	Sample Plan	Name		# of Locations		•	Location
to display the Sample Plan Screen.		,	🕀 Add Ne	w Sample Plan		123	Preset
						=,	Sample Plans
			No Sample	Plans Added		.	Certifications
Add New Sample Plan						*	Users
Press the Add New Sample Plan button to						器	Connectivity
display the New Sample Plan Screen.		Edit		Î	Delete	\$	Instrument
					0 / 50	0	
	♠ 🔳		8			-	Service

Add New Location to a Sample Plan

Press the Add New Location button to display the Sample Plan Setup Dashboard.

New Sample Plan Screen

×	New	/ Sample	Plan		
Group Name	Location	Name	Manual / Preset	2	Location
	Add Nev	w Location		G	Output
				-	Name
	No Da	ata			
Z Edit			📋 Delete		

Sample Plan Setup Dashboard – No group selected


If no group is selected (note the Group Name in this example has a dark grey background color) then all locations will display and be selectable.

If a selected location is linked to a preset

then the Preset button will highlight blue and the manual sample settings button will become disabled. Only the preset linked to the location will appear in the Preset Name list.

× Sample Settings Dashboard						•
Group Name		Location			123	
Group 1	ى	Cocation 0001	\uparrow	-		
Group 2	ى	Cocation 0002		Preset Na	ame	
Group 3	ى	Contention 0003		Preset 1		
	L	ocation 0004				
	L	ocation 0005				
	L	ocation 0006				
	L	ocation 0007				
	L	ocation 0008				
	L	ocation 0009	\checkmark			
3	/ 250		250 / 2500		!	5 / 500

Sample Plan Set Up Dashboard - Group name selected



Sample Plan Setup Dashboard


Select Locations for Sample Plan

Select a Location from the list and it will highlight blue indicating it has been selected.

Manual - Sample Settings

Click on the Manual button and it will highlight blue indicating it has been selected.

Select the Manual Sample Settings button to display the Manual Sample Settings screen. -



Sample Plan Manual Sample Settings

Sample

00:00:00

00:01:00

00:00:00

m³

0.028

Manual

Enter Va

Min Volume: 0.0028

Max Volume: 40.7746

Delay Time


Select the Delay button and enter the delay time using the keypad and press the green Apply button.

Sample Time

Select the Sample Time button and enter the _ sample time using the keypad and press the green Apply button.

Hold Time

Select the Hold button and enter the hold time using the keypad and press the green Apply button.



Sample Time

Mode


Delav


Hold

Cycles Units

Volume

Sample Time





248083466-1 R11

8

Sampl

Chann

Environmen / Ale

Cycles

Select the Cycles button and enter the number of Cycles using the keypad and press the green Apply button.

Units

Select the Units button and select the units, ft3 (cubic feet), m3 (cubic meters), or L (liters).

Select the Volume button and enter the sample volume with the keypad and press the green Apply button.

Number of Cycles



Units of volume



Press the Channel/Alerts buttons to display the Channel Warning and Alarm input screen.

Enable or Disable Channel Sizes

Each particle size channel can be enabled (green) or disabled (grey). If a particle size channel is disabled then no data will be reported on that particle size channel.

Green = Enabled Channel

Grey = Disabled Channel

Note: If Lock Global Channels is enabled then all channels are locked. Additional channels may not be enabled or disabled.

	×	Manual								
	Mode	Sample		Min Samples / Location: 1 Enter Value						Sample
	Delay	00:00:00		Enter va	liue					Channel / Alerts
	Sample Time	00:01:00		1	2	3	delete	ľ		74010
	Hold	00:00:00		4	5	6				
	Cycles	1		,	3					
/	Units	m³		7	8	9				
/	Volume	0.028			0		Apply			

Sample Plan Manual Sample Settings




Manual						
Wan	ning (ft²)	Al	arm (fP)	8 Sample		
	0.0		0.0	Channel Alerte		
	0.0			Environmenta / Alerts		
	0.0		0.0	/ ABETIS		
	0.0					
	0.0		0.0			
	0.0		0.0			

×				
	Channel	Warning (ft ³)	Alarm (ft³)	Sample
	0.3	0.0	0.0	Channel / Alerts
	0.5	0.0	0.0	Environmental / Alerts
	1.0	0.0	0.0	/ Alens
	3.0	0.0	0.0	
	5.0	0.0	0.0	
	10.0	0.0	0.0	

Enter Warning and Alarm Thresholds

Press a Warning or Alarm threshold button to enter particle count threshold values for that particle channel size. Press the Enable button to enable the threshold setting. Or press the Disable button to disable the threshold.



Manual Sam	ple Plan: Cha	nnel, Warning	and Alarms

Red button color indicates the	X Manual					
channel Alarm threshold is enabled.	Channel	Warning (ft ³)	Alarm (ft ³)	8 Sample		
Orange button color indicates the	0.3	0.0	0.0	P Channel / Alerts		
channel Warning threshold is enabled.	0.5	10.0	20.0	Environmental / Alerts		
Dark grey button indicates channel disabled.	1.0	0.0	0.0	/ Alens		
Dark grey button indicates channel disabled.	3.0	0.0	0.0			
Light grey button indicates threshold disabled.	5.0	0.0	0.0			
	10.0	0.0	0.0			

Press the green save button to save the Manual Sample settings to the highlighted location and redisplay the Sample Plan Setup Dashboard.

•

Set Environmental Alerts

Select the Environmental Alerts tab to set up Alarms and Warnings for Environmental Alerts.

Choose Temperature Units

Press to select either Celsius or Fahrenheit. The ApexZ will automatically convert the values within the Alert boxes.

×	< Manual								
Environmental	Units		Alarm High	Warning High	Warning Low	Alarm Low	Sample		
TEMP	°C	۰F	o	0	0	o	Channel Alerts		
RH	q	%	0	0	0	0	Environmenta / Alerts		

Environmental Alert Settings

Set Temperature Alert Values

Press the value for the Alarm High to bring up the input keyboard. Choose Enable or Disable for the Alert, enter a value (within the range indicated) and press the Apply button.

×	5	i5	
Minimum: 0 Maximum: 150	Enable	Disable	
Enter Alarm H	ligh		55
1	2	з	delete
4	5	6	
7	8	9	2
	0	1947). 1	Apply

× 32							
Vinimum: 0 Vaximum: 19	Enablo	Disable					
Enter Warn	ng High		50				
1	2	3	delete				
4	5	6					
7	8	9	at.				
	0	-	Apply				

	32					
	Disable	inable				
			inimum: 0 laximum: 150			
4(Enter Warning Low			
delete	3	2	1			
	6	5	4			
	9	8	7			
Apply		0	1			

×	3	2	
	Enable	Disable	
Minimum: 0 Maximum: 15	o		
Enter Alarm	Low		35
٦	2	3	delete
4	5	6	
7	8	9	
	0	14	Apply

Temperature Alert Settings

Repeat the action for the Warning High, Warning Low, and Alarm Low values.

Set Relative Humidity Alert Values

Press the value for the Alarm High to bring up the input keyboard. Choose enable or Disable, enter a value (between 0 and 100 percent) and press the Apply button.

Repeat the action for the Warning High, Warning Low, and Alarm Low values.

×	< Manual									
Environmental	Units		Alarm High	Warning High	Warning Low	Alarm Low	8	Sample		
TEMP	°C	٩F	55	50	40	35	÷	Channel / Alerts		
RH	c	%	80	70	40	30	8 E	nvironmental / Alerts		

Environmental Alerts Set

×	0465-00553 X	0		×		0		×	Alexand M	0		>>		
	Enable	Disable			Enable	Disable			Enable	Disable				Enable
mum: 0 imum: 100				Minimum: 0 Maximum: 100				Minimum: 0 Maximum: 100		5		Minim Maxim	um: 0 1um: 100	
er Alarm H	9h		80	Enter Warning	High	_	70	Enter Wamin	g Low		40	Ente	r Alarm Low	
1	2	а	delete	١	2	з	delete	1	2	з	delete		1	2
4	5	6		4	5	6		4	5	8			4	5
7	8	9	-	7	8	9		7	8	9			7	8
	ő		Apply				Anniu		0		Apply			0

Relative Humidity Alert Settings

×	()	
	Enable	Disable	
Minimum: 0 Maximum: 10	2		
Enter Narm I	Law		
	2	з	dalete
	5	đ	
7	8	9	
	0	-	Apply


248083466-1 R11

Lighthouse Worldwide Solutions

Add location with settings to the Sample Plan

Press the green Save button to save the highlighted location with settings to the Sample Plan.





Manual Sample Plan Dashboard

Sample Plan with the first location added

In this example "Location 0004" has been added to the new sample plan.

Repeat the above process to add additional locations with Manual or Preset sample settings.

×	New	/ Sample Pl	an	_	
Group Name	Location	Name	Manual / Preset	•	Location
	🕂 Add Nev	v Location			
n/a	Location 0004	Manu	ıal	Ċ	Output
				•	Name
Edit			Delete		

Preset - Sample Settings

Preset sample settings can be assigned to locations in a sample plan.



Select a location and it will

highlight blue indicating that it has been selected.

Select a preset name and it will highlight blue indicating it has been selected.



Press the green Save button to save the selected location with the selected preset sample settings to the new sample plan.

If a location is linked to a preset then the user will only be able to use that preset sample settings at





that location. If a location is not linked to a preset then the user can chose any preset or set up manual sample settings. Note if the location selected already has a Preset assigned to it then only that preset name will be listed on the preset name list. This is because a location may be assigned only 1 preset setting.

Press the preset button and it will highlight blue indicating it has been selected.

Sample Plan with 3 locations added

Example: 3 locations assigned to a Sample Plan.

The first location has manual sample settings and the second two locations have Preset sample settings.

× New Sample Plan								
Group Name	e Location Name	Manual / Preset	•	Location				
	🕀 Add New Loca	ation		_ (12)///				
n/a	Location 0004	Manual	e	Output				
n/a	Location 0005	Preset 1	-	Name				
n/a	Location 0006	Preset 2	•					
	Edit	Delete						

Sample Plan Output Print Options

New Sample Plan

Print

Output


Press the Output button to display the Print Output options:

Enable "On Sample" to print after every sample has been completed.

Enable "Average" to print the Average counts after every sample is completed.

Enable "On Alarm" To print after every alarm.

Press the Name button to display the Sample Plan Name screen.



B

Output

• On Sample

On Alarm

+ Average

Sample Plan Name Screen

Location


Outpu

Name

G

Save the Sample Plan Name

- •
- Press the green Save button to save the new sample plan.



Settings Screen: Sample Plans



The new sample plan name is now listed on the Sample Plans screen.

The total number of locations included in the sample plan is listed in the # of Locations column.

Edit Sample Plan

Select a Sample Plan to edit and it will highlight blue indicating it has been selected.

Press the Edit button to edit the sample plan.


nt	Sample	e Plan N	ame			# of	Locations			•	Location
				÷	Add Ne	w Sample	Plan				
	Sample Plan	3					3			123	Preset
	Sample Plan	2					3			, s	Sample Plans
	Sample Plan	1					3				
										- 1	Certifications
										*	Users
										윪	Connectivity
			Edit				ī	Delete		*	Instrument
									3 / 500		
	^ 1			٠						~	Service

Sample Plans Screen

Edit Sample Plan Locations Screen

Sample Plan 2 ٦ Press the Location button to edit the locations Manual / Preset Group Name Location Name 0 Locatio assigned to the sample plan. Add New Location Location 0004 Manual Group 002 Group 002 Location 0005 Manual Name Select a Location to Edit. Group 002 Location 0006 Manual 🗾 Edit 📋 Delete

Edit Sample Plan Locations Screen



Edit the Manual or Preset Sample Settings

Press the Manual Settings button to select Manual Sample Settings.

Then press the sample setting button to edit the manual sample settings.

Manual Sample settings


Manual -Mode Min Volume: 0.0028 Sample Sample 8 Max Volume: 40.7746 Enter Value Delay 00:00:00 Channel Alert Sample Time 00:01:00 delet Environment: Aler Hold 00:00:00 Cycles Units m³ Volume Apply 0.028

Manual Sample Settings

Manual Channel/Alert settings

Edit Manual Sample Settings.

Manual Warning/Alarm Settings



Warning Threshold


×	(ט	
	Enable	Disable	
Enter Warning			10
1	2	3	delete
4	5	6	
7	8	9	
	0		Apply
			N

Alarm Threshold

×		0	
	Enable	Disable	
Enter Alarm	י		20
1	2	3	delete
4	5	6	
7	8	9	
	0		Apply

248083466-1 R11

Environmental Alerts





Preset Sample Settings

Select Preset for this location.

Press the green Save button to save the Preset sample setting at the selected location.

×	San	nple Plan Setup Das	shbo	ard	•
Group Name		Location		<u>.</u>	123
Group 001		Location 0004		-	
Group 002				Preset Na	ime
Group 003				Preset 1	
				Preset 2	
				Preset 3	
	3/250	1	2500		3 / 500

Select Preset for Location

Edit Sample Plan Output

Press the On Sample button to enable output to the printer after every sample completes.

Press Average to print the average count of the total cycles in a sample.

Press On Alarm to output to the printer after every Alarm threshold is exceeded.

Select Output on Sample or on Alarm

×			Sample	Plan 2		
		в	Print		•	Location
G	Output	•• On Sample			G	Output
		└→ ÷Aver ↓ On Alarm	age		•	Name

Edit Sample Plan Name

Press the Name button to edit the Sample Plan name.

Edit the name using the keypad.

Press the green Apply button to apply the edited Sample Plan name.

Press the green Save button to save the edited Sample Plan settings.

Sample Plan Name Screen




Delete Sample Plan

Select the Sample Plan to be deleted and it will highlight blue indicating it has been selected.

Press the Delete button.

Sample Plan Screen



Press the green check box to confirm deleting the selected Sample Plan.

Or press the red X button to cancel.

Confirm Deleting Sample Plan



Blank Page


Certifications

Add a Certification

From the Home screen press the Configuration button and then press the Certifications button to display the Certifications screen.

Select the + Add New Certification button to display the New Certification input screen.

Settings: Certifications Screen



Certification Types

Select the certification type.

ISO 14644-1: 2015 EU GMP Annex 1 ISO 14644-1: 1999 Federal Standard 209E

New Certification Input Screen

× New Certification						
Туре	ISO 14644-1: 2015		ISO 14644-1: 2015	ŧ	Туре	
Class	5		EU GMP	•	Channel	
Room Status	in-operation				Sample	
Area Units	ft²		ISO 14644-1: 1999	8		
Area	10.000		FED 209E	è	Location	
				-	Name	

Class

Select the cleanroom Class from the classes available for the certification type selected.

ISO 14644-1: 2015 (available classes: 2.0, 2.5, 3.0, 3.5, 4.0,


4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9)

×	New Ce	rtif	ication				
Туре	ISO 14644-1: 2015		2	2.5	个	Ħ	Туре
Class	5	İ	3	3.5		••	Channel
Room Status	in-operation	i	4	4.5		8	Sample
Area Units	ft²					•	Location
Area	10.000	1	5	5.5		Ó	
		I	6	6.5			Name
			7	7.5			
					\mathbf{A}		

Cleanroom Class


Certification Classes & Grades

The available classes/grades will change according to the certification type selected above.



Room Status


Certifications will require Room Status to be selected. Select in-operation, at-rest, or as-built and the selected room status will highlight blue indicating it has been selected.



ISO 14644-1: 2015

Federal Standard 209E

Туре	FED 209E	in-operation	🗮 Тур
Class	100	at-rest	🔶 Channe
Airflow	non-unidirectional		
Room Status	in-operation	as-built	Sampl
Area Units	ft²		Locatio
Area	10.000		Nam



ISO 14644-1: 1999

EU GMP Annex 1

Туре	EU GMP	in-operation	i=	Тур
Grade	А	at-rest		Channe
Room Status	in-operation			
Area Units	ft²		8	Sampl
Area	10.000		2	Locatio
				Nam

EU GMP certifications will be limited to Room Status of In-Operation or At-Rest only.

Airflow

For Federal Standard 209E Certifications, airflow will be required to select either "non-unidirectional" or "unidirectional".

Туре	FED 209E	non-unidirectional
Class	100	unidirectional
Airflow	non-unidirectional	
Room Status	in-operation	
Area Units	ft²	
Area	10.000	

Area Units

New Certification Press the Area Units button. × Туре ISO 14644-1: 2015 ≣ Туре Class 5 Channel Room Status in-operation Select the ft² or m² area units. 8 Sample Area Units ft² Location 9 Area 10.000 Name •

Area

Area

Area Units

Click on the Area button and the numeric	×	New Ce	erti	ficat	ion				
keypad will pop up.	Туре	ISO 14644-1: 2015			Max A		ea: 1.000	≔	Туре
$\langle \rangle$	Class	5	i I	Enter V			·	<u>.</u>	Channel
\sim	Room Status	in-operation	ļ		_				(
	Area Units	ft²		1	2	3	delete	8	Sample
	Area	10.000		4	5	6		è	Location
Enter the Area value.				7	8	9		•	Name
Press the green Apply button.					0		Apply		

Press the Channel button to display the New Certification Channels screen next.

Channels

Certification channel sizes are limited and locked according to the certification type.




Enabled channels will appear with a green check box.



Disabled channels appear with a red x box.

Press the green check box to disable a channel and it will change to a red X box and vise-versa.

Certification Channels



Note the on screen message "Some allowed channels are not enabled".

In this example the 0.3 μ m size particle channel is an allowed channel but it was disabled.

The dark grey null set button indicates this channel size is not included in the certification.

Press the Sample button to display the certification sample settings next.

Sample settings for new Certification.

Delay Time


Select the Delay button and enter the delay time using the keypad and press the green Apply button.


Sample Time

Select the Sample Time button and enter the sample time using the keypad and press the green Apply button.

Hold Time

Select the Hold button and enter the hold time using the keypad and press the green Apply button.





Sample Time

Min Time: 00:10:19 Max Time: 23:59:59 Enter Value							
	00:00:00						
1	2	3	delete				
4	5	6					
7	8	9					
	0		Apply				

Hold Time							
Max Hold: 23:59:59 Enter Value							
	00:00:00						
1	2	3	delete				
4	5	6					
7	8	9					
	0		Apply				

Cycles

Select the Cycles button and enter the number of Cycles using the keypad and press the green Apply button.

Volume

Select the Volume button and enter the sample volume with the keypad and press the green Apply button.



	Volume							
Min Volume: 10.3167 Max Volume: 1,439.9833 Enter Value								
1	2	3	delete					
4	5	6						
7	8	9						
	0	•	Apply					

Select the Location button next to display the certification location selection screen.


Locations

If a Group Name is selected then only those locations assigned to that group will display in the Location List.


Certification: Locations

e locations cation List.	× New Certification									
	🤆 Group		(Group 1	∷	Туре				
	Locations : 3			selected : 0						
	Location Name		Assigned Lo	cations	•	Channel				
	Add Available Location Select All	s			8	Sample				
	Location 0001		No Locatio	n colocted	2	Location				
	Location 0002		No Locatio	IT Selected		Name				
*	Location 0003					, iaine				


If no Group Name is selected the all locations will be available to be assigned to the new certification.



Certification: Locations



Certification Locations



Select one or more location names and each will highlight blue indicating it has been selected. Each location selected will be displayed under the Assigned Locations column.

To undo a selected location from the assigned list simply press the location name to de-select and it will highlight grey and that location will be removed from the assigned locations list.


Press the Name button to display the Enter Certification Name screen next.

Enter the New Certification name using the keypad. Press the green Apply button to apply the new certification name.



Certification Name - Apply

Certification Name




Settings: Certifications Screen




Edit a Certification

From the Home screen press the Settings button and then press the Certifications button.



Edit certification settings


Select a certification setting button to edit for:



Delete a Certification


selected.

From the Home screen press the Settings button and then press the Certifications button.



Confirm Delete Certification

Confirm Configuration Delete



Blank Page

Users & Permission Levels

Users

ApexZ supports up to 1000 users and 10 user Permission Levels. Username and Password can be kept locally on the ApexZ unit or verified through Active Directory on your company network. ApexZ supports assigning Permission Levels to individual Active Directory Groups. In addition ApexZ can prioritize AD Groups and their respective permission levels.

Permission Level User Name • Loca Add New Permission Levels Here -Add New Permission Level Add New User 123 Pres Administrator Add New Local Users Here Power User Sample Plans Operator No Users Added Certification Unassigned ---Users Connectivity 🗾 Edit Delete 📝 Edit Delete Ż. Instrumen 4/10 3 / 1000 -Servi

User Name and Permission Levels Screen

ApexZ supports up to 10 Permission Levels.

Administrator Power User Operator Unassigned (only for unassigned users)

... up to 10 total permission levels.

User permissions can be enabled or disabled for each permission level created.


Operator and Power User Permission Levels can be edited or deleted.

The "Administrator" and "Unassigned" permission levels may not be edited or deleted.

Permission Levels Settings

×	Operator 5										
Main Sc	reen	Data Reco	ord Screen	Data Repo	rt Screen	Setting 5	icreen	🔩 Pe	rmission	s	
Start / St	ор						 Image: A second s	<u></u> A	D Group	s	
🔲 Data Rec	ord Scree	'n					× .				
🖹 Data Rep	ort Scree	n					× .	C≣	Priorit	У	
🔅 Setting 🕈	×			New	Permiss	sion Lev	el				
	Mair	n Screen	Data Rec	ord Screen Data Report Screen			Settin	g Screen	**	Permission	
	Filter	Data Record			<u> </u>	1	Users				
	Delete All Data Records							~	-		
	Sec. 1	💊 Add Tag							¢≡	Priority	
	🖹 Add I	Report						~		Name	
	G*Export Data Records								•		
	Ь ∲ ∪ѕв							1			
	L}							×			
								V			
		> 🗟 Networ	k Printer					V .	\mathbf{r}		


Add a New User



User Type

User Type

User type "local" will verify Username and password saved locally on the ApexZ instrument.



Note: When creating an Active Directory User no password is required. The password resides on the Network Active Directory.

Note: To enable Users, at least one Admin level user must be created first.

Username

Enter user name on keyboard.


Press the green apply button and the new username will display in the username box.



Next Press the Password button.



Username Screen



Password Screen

Password

Enter Password on the keyboard.

New User × Press the green apply button to confirm password entered. Туре Username Password Password Re-enter the same password to Name confirm password then press the green apply button again to confirm. k Permission Level Shift Delete Space Next Press the Name button.



Click here to toggle view or hide password.

First and Last Name

Press the Name tab below the Password button to advance to the First and Last Name screen.

Enter First Name from the keyboard.

Press Apply and the first name will be displayed above.

Enter the Last Name from the keyboard.

Press Apply and the last name entered will be displayed in the last name box.

After successfully applying the first and last name; Press the Permission Level button to select the permission level.

Permission Level

Select a Permission Level; Administrator, Power User or Operator. The Permission Level will highlight blue indicating it has been selected.

Press the green save button to save the Permission Level and display the new user included on the user screen.

×		Operator 1											
[Peter]					Smith			Туре
First Na	ame	me Peter					Last Name Smith			Smith	2	Username	
	q	w	е	r	t	у	ŭ	i	o	р		٩.	Password
		а	s	d	f	g I	h j	k	Т			•	Name
	Shift	z	x	с	v	b	n	m		Delete		*	Permission Level
	123				Sp	ace			Ap	ply			

First and Last Name Screen

Select Permission Level Screen


×		Operator 1		•
	#	Permission Level		Туре
	1	Administrator	~	Username
	2	Power User		
	3	Operator	٩	Password
	4	Unassigned		Name
			~ 2	Permission Level

Edit User


To edit users, an Administrator level user must be logged in or Users must be disabled.

From the User screen select a User to edit from the list and the user name will highlight blue to indicate it has been selected.

Press the Edit button to display the User setting tabs.



Users Screen




Note: Administrator level users can edit their own password but cannot edit other user passwords. If other users forget their password then a logged in Administrator Permission Level User will need to recreate their account.

248083466-1 R11

Users Screen

Delete a User

Users Screen



User Login


When Users are enabled ApexZ will require a Username and Password logon.

Enter Username and Password


Touch the Username box and the border will highlight yellow indicating it has been selected.

If Users are enabled then login with your local ApexZ Username and password.

If Active Directory is enabled then login with your network Active Directory username and password.



Login Screen



Home Screen

Admin 1 Ų۶ 06/18/2019 2:24:27 pm 04 X 0 Group 1 0 Location 0001 ¢ Σ μm 0 0.3 Sample Plans 0.5 0 1.0 0 Manual Sample 3.0 0 :00:00 0 5.0 00:00:00 9 / 10000 10.0 0 IDLE 📑 Feed

The username will display in the username box as you type.

Touch the Password button and the border will highlight yellow indicating it has been selected.

Enter your user password with the keypad.


The login button will illuminate green indicating it is ready to be selected.

Press the green Login button to login to the ApexZ.


The logged in Username will be displayed on the Home page.

ApexZ Idle


User Logout



Logout Screen when ApexZ is idle




Logout Screen when ApexZ is sampling



Administrator Permission Level

Administrator Permission level users have all ApexZ features enabled and no features may be disabled.

Add new Permission Level



User Screen – Add New Permission Level

Edit "Power User" Permissions Screen


Select the screen tabs to edit permission levels. • × Power[']User Main Screen Data Record Screen 喿 Permissions Main Screen Data Record Screen Data Report Screen Setting Screen Data Report Screen Setting Screen Start / Stop Users * Data Record Screen Priority c= 🗕 Data Report Screen 🙁 Setting Screen Name


Edit Permission Levels


Select the screen tab to edit permissions for that screen. Enable or disable any feature by pressing the button.



Press the green Save button to save the Permission Level settings.







Assign AD Group to Permission Level

Active Directory Group names can be manually assigned to a Permission Level or by logging into the network and selecting from a list of AD Groups.

Manually Assign Permission Level Press the + Manual button to display the Enter AD Group Name screen. AD Groups Assigned AD Gr

Assign Active Directory Group to Permission Level

Enter Active Directory Group Name

Enter the AD Group Name to assign to the permission level.

Press the green Apply button.

D Group	Name	40-		44					+2
q	w	е	r	t	у	u	i		р
	а	s	d	f	g h	j	k	I	
Shif	z	×	с	v	ь	n	m		Delete
123				Sp	ace			Ap	ply

Login to select AD Groups and assign Permission Levels

Note: To login and display network AD Groups the AD URL and Domain must be entered in the ApexZ Settings Connectivity AD/LDAP screen.

Press the Login button to login to the Active Directory.



Assign Active Directory Group to Permission Level

Enter the Active Directory Username and Password.




Press the Apply button.

After a successful login AD Groups will display in the AD Groups list column.

Active Directory: Username and Password




Assign 1 AD Group to 1 Permission Level




Press an AD Group name and it will highlight blue and display in the "Assigned AD Groups" column.

This AD Group has now been assigned to this Permission Level.



Prioritize Permission Levels

ApexZ can prioritize the permission levels created. This is because some users may appear in more than one Active Directory Group. If a user login appears in more than one Active Directory group (with different permission level settings) then the permission level assigned will be from the higher priority permission level setting on the Priority Screen.



Default Permission Levels


Permission levels:	Administrator	Power User	Operator	
	(Enable Locked)	(Enable/Disable)	(Enable/Disable)	✓ Enabled
Main Screen				
Start & Stop	✓	\checkmark	\checkmark	
Data Screen	\checkmark	\checkmark	\checkmark	
Report Screen	\checkmark	\checkmark	\checkmark	
Settings Screen	✓	\checkmark		
Data Screen	Administrator	Power User	Operator	Disabled
Filter Data Records	\checkmark	\checkmark	\checkmark	
Delete All Data Records	\checkmark	\checkmark	•	
Add Tag	\checkmark	\checkmark	\checkmark	
Add Report	\checkmark	\checkmark	\checkmark	
Export Data Records	\checkmark	\checkmark	\checkmark	
Export – USB	\checkmark	\checkmark	\checkmark	
Export – Print	\checkmark	\checkmark	\checkmark	
Export – CIFS	\checkmark	\checkmark	\checkmark	
Export – Network Printer	\checkmark	\checkmark	\checkmark	
Data Report Screen	Administrator	Power User	Operator	
Filter Data Reports	\checkmark	\checkmark	\checkmark	
Delete All Data Reports	\checkmark	\checkmark		
Delete Data Report	\checkmark	\checkmark		
Export Data Reports	\checkmark	\checkmark	\checkmark	
Export Report – USB	\checkmark	\checkmark	\checkmark	
Export Report – Print	\checkmark	\checkmark	\checkmark	
Export Report – CIFS	\checkmark	\checkmark	\checkmark	
Export Report – Network Printer	\checkmark	\checkmark	\checkmark	
Filter Audit Trails	\checkmark	\checkmark	\checkmark	
Export Audit Trails	\checkmark	\checkmark	\checkmark	
Export Audit Trails – USB	\checkmark	\checkmark	\checkmark	
Export Audit Trails – CIFS	\checkmark	\checkmark	\checkmark	
Set Audit Trail Comments	\checkmark	\checkmark	\checkmark	
Setting Screen	Administrator	Power User	Operator	
Sample Parameters	\checkmark	\checkmark		
Users	\checkmark			
Ethernet / Wi-Fi	\checkmark	\checkmark		
AD/LDAP/CIFS/Network Printer/	\checkmark	\checkmark		
3 rd Party API				
Time	\checkmark	\checkmark		
Audio/Visual	\checkmark	\checkmark		1
Options	\checkmark	\checkmark		1
Output	\checkmark	\checkmark		
Configuration	✓	√		
Service	\checkmark	\checkmark]

Service

From the Settings screen press the Service button to view the Maintenance, About and Service Logon tabs.

Maintenance

Press the Maintenance button to display the Firmware, Reset and Calibration functions.



Press the green check button to confirm and update the firmware on the ApexZ from the USB key.

Or press the red X button to cancel.



Instrument Settings Reset

Press the Instrument Settings Reset button to reset all ApexZ instrument settings.

Locatio About 9 Maintenance Login Firmware Pres Update 🕒 Update Firmware Sample Plan Reset Certification Reset Instrument Settings Reset Factory Reset Calibration Connectivity 1 Month set Calibration Reminder \$ Instru . Servic

Service: Instrument Settings Reset


Settings Reset

Or press the red X button to cancel.

Factory Reset

Pressing Factory Reset will DELETE all data on the ApexZ and return the instrument's configuration to its original factory default configuration. Users will be reset to the default "Administrator", "Power User" and "Operator" Permission Levels. There is no data recovery after a Factory Reset. All user created data will be deleted.

Press the Factory Reset button to prepare to reset to factory settings.




Locatio Maintenance About Login 123 Firmware Prese Update Update Firmware Sample Plar Reset 👯 Certification Reset Instrument Settings Reset Sactory Reset -Calibration Connectivity 1 Month set Calibration Reminder ¢. Instrum . Servic

Service: Factory Reset

Press the green check button to confirm factory reset or press the red X button to cancel and NOT factory reset.

Calibration Reminder

Press the Calibration Reminder button to select the time period for a reminder status warning before the calibration due date. An orange warning will appear in the Self-Diagnostices Status display.



About Screen

Instrument Name

Press the Instrument Name button to enter a name for your ApexZ instrument.


Use the elevator bar to scroll down and view the Instrument identification.

Model Serial # Sensor # Calibration Due Date

View the ApexZ firmware.

MEAN Firmware Version MAIN Firmware Version SIB Firmware Version POWER Firmware Version

Firmware



About Instrument

About: Firmware & Open Source Libraries

Maintenance	About	Login	•	Location
Cal Due:			123	Preset
	Firmware			
JS MEAN Firmware		2.02.010	_,	Sample Plans
C MAIN Firmware		0.01.001	EM	Certifications
SIB Firmware		1.01.003		
UPOWER Firmware		1.01.004	-	Users
Open	Source Libraries		品	Connectivity
♦ Open Source Libraries		Libraries	\mathbf{v}	instrument
♠ ≡ 🖹 🌣			٩	Service

Open Source Libraries

Open Source Libraries

Press on the Libraries button to view the open source libraries used by ApexZ.

Licenses

Libraries We Use	
Thank you to the open source community!	
Node-ActiveDirectory	\uparrow
Ag-Grid	
Another JSON Schema Validator (Ajv)	
Advanced Linux Sound Architecture (ALSA)	
Advanced Linux Sound Architecture (ALSA) utilities	
AngularJS	
AngularJS Bootstrap	
AngularJS ChartJS	
AngularJS Local Storage	

Service Login

Press the Settings then Service button to display the Service Login screen.

Note the ApexZ Passcode number and use this code to get a password from Lighthouse Service Department calibration website.

Service Login Use the keypad and enter the password from the LWS Service Department here. Maintenance About Login 030495 Retries Passcode: 123 Pre Press the green Apply button. Enter Password: 🛼 Sample Plar Follow ApexZ Service Manual for instructions on how Certificat to change ApexZ service parameters. ----User 윦. Connectivit Space \$.

Calibration

Calibration – Contact Authorized Facility

ApexZ3 is factory calibrated to 1 CFM (1 cubic foot per minute).

ApexZ50 is factory calibrated to 100 LPM (100 liters per minute)(3.53 cubic feet per minute).

Contact a factory authorized facility to calibrate the ApexZ.

Calibration – Due Date

To display calibration due date, see Section: Instruments / About.

ApexZ calibration is valid for 1 year from date of calibration.

After the ApexZ calibration data has passed, all records will be shown with a strike-through on the data to indicate that the data is not valid.

Trouble Shooting

Printer:

Printer paper not advancing while printing could be the printer door is not closed.

- If the printer door is not closed then the paper will not advance but the printer will continue to print. Quickly close the printer door and printing will advance the paper properly.

Printer printing but no ink is printing on the paper could be the printer paper is loaded backwards.

If the printer paper is load backwards the paper may still advance forward but no report data will
print on the paper. After printing has stopped, simply open the printer door, remove the paper roll
and reverse so that the paper feeds up from the back side of the roll and out the top of the printer.
(see section: Printer Paper Installation)

Printer paper looped back into the printer housing!

- If the advancing printer paper loops back into the printer housing it may cause a paper jam.
- If this happens it is recommended to stop the printing. Simply open the printer door at any time and printing will pause. You can adjust the printer paper and close the printer door.
- Printing will resume where it left off.

Printer paper not straight

- If the printer paper is not straight or misaligned simply open the printer door and printing will stop.
- Realign the printer paper so it is feeding straight and correctly.
- Line up the left and right edges with the printer door exit slot.
- Close the printer door. If printing was interrupted then it will resume.
- See Section: Printer Paper Installation

Clean and Disinfect

External Enclosure

The following cleaning agents are acceptable for cleaning ApexZ polycarbonate exterior surfaces.

Spor Klenz
Bleach
Hydrogen Peroxide
IPA

Warnings

To minimize the chance of electric shock, turn off and disconnect the power adapter during sterilization. Do not submerge ApexZ in any liquids.

Do not spray disinfectant solution directly into ApexZ enclosure.

Do not use gases to disinfect ApexZ.

Lighthouse will VOID the warranty if an unsuitable disinfectant solution damages the ApexZ instrument.

Terminology

Button:	Area on the touch screen the user can press and ApexZ will accept input and commands.
Enable:	Green box with yellow check mark: this button indicates enabled.
Disable:	Light gray box: this button indicates disabled.
Enter with the keypad:	Touch the screen keypad with your finger and type the required input. You may need to press the green Apply button after entering text or numbers.
Inactive Button:	Button has a dark background color indicating the button is inactive and will not respond when pressed.
Press:	Touch the screen with your finger.
Select:	Touch the screen with your finger.
Toggle button:	This is a button that displays only 2 options, touch or slide to change between options.

Appendix A: Limited Warranty

Limitation of Warranties

A. Lighthouse Worldwide Solutions (LWS) warrants that all equipment shall be free from defects in material and workmanship under normal use for a period of two years from date of shipment to Buyer. LWS does not warrant that operation of the software will be completely uninterrupted or error free or that all program errors will be corrected. Buyer shall be responsible for determining that the equipment is suitable for Buyer's use and that such use complies with any applicable local, state, or federal law. Provided that Buyer notifies LWS in writing of any claimed defect in the equipment immediately upon discovery and any such equipment is returned to the original shipping point, transportation charges prepaid, within two years from date of shipment to Buyer and upon examination LWS determines to its satisfaction that such equipment is defective in material or workmanship, i.e. contains a defect arising out of the manufacture of the equipment and not a defect caused by other circumstances, including, but not limited to accident, misuse, unforeseeable use, neglect, alteration, improper installation, improper adjustment, improper repair, or improper testing, LWS shall, at its option, repair or replace the equipment, shipment to Buyer prepaid. LWS shall have reasonable time to make such repairs or to replace such equipment. Any repair or replacement of equipment shall not extend the period of warranty. If the Instrument is modified or in any way altered without the explicit written consent of LWS then the warranty is null and void. This warranty is limited to a period of two years, except as noted below, without regard to whether any claimed defects were discoverable or latent on the date of shipment. The length of warranty for pumps in hand held particle counters is one (1) year. Batteries and accessories with all products are warranted for one (1) year. Fuses and purge filters carry no warranty. If a third party battery is used in the product, the product warranty is null and void. If the battery is charged by a third party battery charger the battery warranty is null and void.

B. If Buyer shall fail to pay when due any portion of the purchase price or any other payment required from Buyer to LWS under this contract or otherwise, all warranties and remedies granted under this Section may, at LWS's option, be terminated.

C. THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER REPRESENTATIONS, WARRANTIES AND COVENANTS, EXPRESS OR IMPLIED WITHRESPECT TO THE EQUIPMENT AND ANY DEFECTS THEREIN OF ANY NATURE WHATEVER, INCLUDING AND WITHOUT LIMITATION WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. LWS SHALL NOT BE LIABLE FOR, AND BUYER ASSUMES ALL RISK OF, ANY ADVICE OR FAILURE TO PROVIDE ADVICE BY LWS TO BUYER REGARDING THE EQUIPMENT OR BUYERS USE OF THE SAME. UNDER NO CIRCUMSTANCES SHALL LWS BE LIABLE TO BUYER UNDER ANY TORT, NEGLIGENCE, STRICT LIABILITY, OR PRODUCT LIABILITY CLAIM AND BUYER AGREES TO WAIVE SUCH CLAIMS. LWS'S SOLE AND EXCLUSIVE LIABILITY AND BUYERS SOLE AND EXCLUSIVE REMEDY, FOR ANY ON CONFORMITY OR DEFECT IN THE PRODUCTS OR ANYTHING DONE IN CONNECTION WITH THIS CONTRACT, IN TORT, (INCLUDING NEGLIGENCE), CONTRACT, OR OTHERWISE, SHALL BE AS SET FORTH IN THE SUBSECTION A HEREOF AS LIMITED BY SUBSECTION B HEREOF. THIS EXCLUSIVE REMEDY SHALL NOT HAVE FAILED OF ITS ESSENTIAL PURPOSE (AS THAT TERM IS USED IN THE UNIFORM COMMERCIAL CODE) PROVIDED THAT THE SELLER REMAINS WILLING TO REPAIR OR REPLACE DEFECTIVE EQUIPMENT (AS DEFINED IN SUBSECTION A) WITH A COMMERCIALLY REASONABLE TIME AFTER RECEIVING SUCH EQUIPMENT. BUYER SPECIFICALLY ACKNOWLEDGES THAT SELLER'S PRICE FOR THE EQUIPMENT IS BASED UPON THE LIMITATIONS OF LWS'S LIABILITY AS SET FORTH IN THIS CONTRACT.

Warranty of Repairs After Initial Two (2) Year Warranty

A. Upon expiration of the initial two-year warranty, all parts and repairs completed by an authorized Lighthouse repair technician are subject to a six (6) month warranty.

B. Other than the above, LWS makes no warranty of any kind expressed or implied, except that the products manufactured and sold by LWS shall be free from defects in materials and workmanship and shall conform to LWS's specifications; Buyer assumes all risk and liability resulting from use of the products whether used singly or in combination with other products. If instrument is modified or in any way altered without the explicit written consent of LWS, then the warranty is null and void.

C. WARRANTY REPAIRS SHALL BE COMPLETED AT THE FACTORY, BY AN AUTHORIZED SERVICE LOCATION, BY AN AUTHORIZED SERVICE TECHNICIAN, OR ON SITE AT BUYER'S FACILITY BY A LIGHTHOUSE AUTHORIZED EMPLOYEE. BUYER PAYS FREIGHT TO FACTORY; SELLER WILL PAY STANDARD RETURN FREIGHT DURING THE WARRANTY PERIOD. BUYER MAY SELECT A FASTER METHOD OF SHIPMENT AT ITS OWN EXPENSE. This Page Intentionally Left Blank.

Service and Support Tel. 1-800-945-5905 (USA Toll Free) Tel. 1-541-770-5905 (Outside of USA) techsupport@golighthouse.com www.golighthouse.com

