

# SHIELD-P2-PRO

## FILTERING HALF MASK



Leading the world in contamination monitoring has inspired us to create the most innovative FFP2 non-reusable filtering face piece on the market. The three-panel design naturally contours to your face to keep the soft nose foam seal in place.

The SHIELD-P2-PRO is a product of countless iterations and focuses on what matters most to the end-user.


### FEATURES

- 92% Fit Test Passing on a Wide Range of Faces
- 25% Increase in Visibility
- Built in Anti-Fogging Technology

### SHIELD-P2-PRO The Fit Test Champion

<b>(EU) 2016/425 Type-Examination Certificate</b>	NIOSH # TC 84A-9466 (EN 149:2001 + A1:2009 FFP2 NR)
<b>AS/NZS 1716:2012 Class P2 Certificate</b>	230889-A(i)
<b>S. Africa Approval</b>	NRCS/8072/0407
<b>Material Composition</b>	
<b>Filter</b>	Polypropylene
<b>Nose Foam</b>	Medical Grade PVC
<b>Nose Clip</b>	Dual Iron Core Polypropylene Coated Wire
<b>Cover web</b>	Polypropylene
<b>Head Straps</b>	Polyisoprene (Black)
<b>Shelf Life</b>	3 Years from the manufacturing date
<b>Storage Conditions</b>	Between -4°F and 86°F Not Exceeding 80% RH
<b>Country of Origin</b>	Manufactured in USA from internationally sourced materials
<b>Manufactured By</b>	Lighthouse Worldwide Solutions



 Request Samples

#### Time Use Limitation:

The respirator may be used until damaged, breathing becomes difficult or contaminated with fluids. Discard after every use. Follow national, state, local, and facility infection control guidance and policies.

#### Storage Instructions:

Store respirators in original packaging and away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals.



**MANUFACTURED IN USA  
FROM INTERNATIONALLY  
SOURCED MATERIALS**



Lighthouse Worldwide Solutions reserves the right to change specifications without notice.

## Who We Are

For more than 40 years, Lighthouse Worldwide Solutions has met the indoor contamination monitoring and control needs of the FDA-regulated healthcare, pharmaceutical, and medical device manufacturing industries. We serve customers in more than 70 nations from our USA-based HQ and multiple international offices. Our customers have leveraged our regulatory and contamination control expertise to help ensure their success.

**Our expertise in contamination control has inspired us to create the most innovative FFP2 disposable respirator to enter the market.**

Lighthouse Worldwide Solutions has developed the SHIELD-P2-PRO disposable respirator to fit more comfortably and securely than any other respirator on the market. The SHIELD focuses on what matters most to the end-user.

Our customers have told us the SHIELD is a “fit test champion,” checking the boxes in fit, comfort, and breathability. In addition to a good seal, we know visibility is also crucial; our patent-pending ClearView Bar can increase your peripheral viewing angle up to 25%. As a NIOSH N95 certified respirator, the SHIELD is the ideal choice for first responders and all others who require a superior quality and high performance FFP2 disposable respirator.

**Our products are made in the USA.**

**We look forward to discussing how we can be of service to your team!**

