

# SHIELD-95-MEDI

S U R G I C A L   N 9 5   D I S P O S A B L E   R E S P I R A T O R



**L** eading the world in contamination monitoring has inspired us to create the most innovative Surgical N95 Disposable Respirator on the market. The three-panel design naturally contours to your face to keep the soft nose foam seal in place.

The SHIELD-95-MEDI is a product of countless iterations and focuses on what matters most to the end-user.

## F E A T U R E S

- 92% Fit Test Passing on a Wide Range of Faces
- 25% Increase in Visibility
- Built in Anti-Fogging Technology

### SHIELD-95-MEDI The Fit Test Champion

<b>CERTIFICATION BY</b>	NIOSH # TC 84A-9431
<b>Material Composition</b>	
<b>Filter</b>	Polypropylene
<b>Nose Foam</b>	Medical Grade PVC
<b>Nose Clip</b>	Dual Iron Core Polypropylene Coated Wire
<b>Cover web</b>	Polypropylene
<b>Head Straps</b>	Polyisoprene (Blue)
<b>Fluid Resistance</b>	ASTM 2100 Level 3 (160 mmHg)
<b>Flammability Rating</b>	Class 1
<b>Shelf Life</b>	3 Years from the Date of Manufacturing
<b>Storage Conditions</b>	Between -4°F and 86°F Not Exceeding 80% RH
<b>Country of Origin</b>	Manufactured in USA from Internationally Sourced Materials
<b>Manufactured By</b>	Lighthouse Worldwide Solutions


#### Time Use Limitation:

The respirator may be used until damaged, breathing becomes difficult or contaminated with fluids. Discard after every use. Follow national, state, local, and facility infection control guidance and policies.

#### Storage Instructions:

Store respirators in original packaging and away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals.



 Learn More



**MANUFACTURED IN USA  
FROM INTERNATIONALLY  
SOURCED MATERIALS**

Lighthouse Worldwide Solutions reserves the right to change specifications without notice.

## Who We Are

For more than 40 years, Lighthouse Worldwide Solutions has met the indoor contamination monitoring and control needs of the FDA-regulated healthcare, pharmaceutical, and medical device manufacturing industries. We serve customers in more than 70 nations from our USA-based HQ and multiple international offices. Our customers have leveraged our regulatory and contamination control expertise to help ensure their success.

**Our expertise in contamination control has inspired us to create the most innovative surgical N95 disposable respirator to enter the market.**

Lighthouse Worldwide Solutions has developed the SHIELD-95-MEDI surgical N95 disposable respirator to fit more comfortably and securely than any other respirator on the market. The SHIELD-95-MEDI focuses on what matters most to the end-user: fit, breathability, visibility, and durability.

Our customers have told us the SHIELD-95-MEDI is a “fit test champion,” checking the boxes in fit, comfort, and breathability. In addition to a good seal, we know visibility is also crucial; our patent-pending ClearView Bar can increase your peripheral viewing angle up to 25%. As an ASTM Level 3 surgically rated / NIOSH N95 certified respirator, the SHIELD-95-MEDI is the ideal choice for healthcare workers and all others who require a superior quality and high performance surgical N95 disposable respirator.

**Our products are made right here in the USA!**

**We look forward to discussing how we can be of service to your team!**

