

# SHIELD-95-PRO

## N95 DISPOSABLE RESPIRATOR



Leading the world in contamination monitoring has inspired us to create the most innovative N95 Disposable Respirator on the market. The three-panel design naturally contours to your face to keep the soft nose foam seal in place.

The SHIELD-95-PRO is a product of countless iterations and focuses on what matters most to the end-user.

### FEATURES

- 92% Fit Test Passing on a Wide Range of Faces
- 25% Increase in Visibility
- Built in Anti-Fogging Technology

### SHIELD-95-PRO The Fit Test Champion

<b>CERTIFICATION BY</b>	NIOSH # TC 84A-9466
<b>Material Composition</b>	
Filter	Polypropylene
Nose Foam	Medical Grade PVC
Nose Clip	Dual Iron Core Polypropylene Coated Wire
Cover web	Polypropylene
Head Straps	Polyisoprene (Black)
<b>Style</b>	3 Panel, Flat-fold
<b>Regulation</b>	Meets NIOSH 42 CFR 84
<b>Shelf Life</b>	3 Years from the Date of Manufacture
<b>Storage Conditions</b>	Between -4°F and 86°F Not Exceeding 80% RH
<b>Country of Origin</b>	Manufactured in USA from Internationally Sourced Materials
<b>Manufactured By</b>	Lighthouse Worldwide Solutions



Request Samples

#### Time Use Limitation:

The respirator may be used until damaged, breathing becomes difficult or contaminated with fluids. Discard after every use. Follow national, state, local, and facility infection control guidance and policies.

#### Storage Instructions:

Store respirators in original packaging and away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals.



MANUFACTURED IN USA  
FROM INTERNATIONALLY  
SOURCED MATERIALS

Lighthouse Worldwide Solutions reserves the right to change specifications without notice.

## Who We Are

For more than 40 years, Lighthouse Worldwide Solutions has met the indoor contamination monitoring and control needs of the FDA-regulated healthcare, pharmaceutical, and medical device manufacturing industries. We serve customers in more than 70 nations from our USA-based HQ and multiple international offices. Our customers have leveraged our regulatory and contamination control expertise to help ensure their success.

**Our expertise in contamination control has inspired us to create the most innovative N95 disposable respirator to enter the market.**

Lighthouse Worldwide Solutions has developed the SHIELD N95 disposable respirator to fit more comfortably and securely than any other respirator on the market. The SHIELD focuses on what matters most to the end-user.

Our customers have told us the SHIELD is a “fit test champion,” checking the boxes in fit, comfort, and breathability. In addition to a good seal, we know visibility is also crucial; our patent-pending ClearView Bar can increase your peripheral viewing angle up to 25%. As a NIOSH N95 certified respirator, the SHIELD is the ideal choice for first responders and all others who require a superior quality and high performance N95 disposable respirator.

**Our products are made right here in the USA!**

**We look forward to discussing how we can be of service to your team!**

